

reed
by eqwal.

Off the shelf

SPRING 2026



EDITION 5

WELCOME

Reed Medical's 'off-the-shelf' footwear range provides ready-made, medical-grade shoes and acute phase off-loading designed to offer protection, comfort and support, and for patients living with painful foot conditions.

The C/E marked medical footwear is suitable for individuals with diabetes, swollen feet, sensitive skin, mild deformities and post-surgery off-loading.

- Diabetic protected activity footwear
- Protective shoes for primary prevention
- Therapeutic shoes for secondary prevention
- Acute phase off-loading

Key features typically include extra depth and width fittings, soft safeseam linings to reduce friction, cushioned insoles for pressure relief, and supportive yet lightweight flexible, semi-rigid, rigid biomechanical soles to aid stability and mobility. Adjustable fastenings also allow for easy fitting and accommodation of swelling or orthotics.

These shoes offer a practical, cost-effective alternative to fully bespoke footwear while still meeting clinical needs for everyday wear.



CONTENTS

Podartis	
Deambulo Line	4-5
Activity Line	6-14
Diabetic & Rheumatic	15-17
Post-Surgery & Acute Phase Off-Loading	18-19
Podowell	20-22
Lymphoedema Shoes	23
Sizing Guide	24-27

DEAMBULO LINE

Deambulo therapeutic shoes have been designed to increase the wearer's stability, prevent falls and accommodate any deformities.



DEAMBULO L

Style no. SA200

Size 39-47

Fitting 12 - Variable volume

Size chart A

- Stretchable upper in Setaflex
- No seams at friction points
- Lightweight anti-stumble flexible outsole
- Extended stabilising back counter



DEAMBULO X

Style no. SA400

Size 36-47

Fitting Variable volume

Size chart A

- Fully opening with variable volume
- Low stabilising back counter
- Biomechanical anti-stumbling outsole
- No seams at friction points
- Super lightweight



DEAMBULO SPORT

Style no. SA300

Size 35-41

Fitting 12 - Variable volume

Size chart B

- Stretchable upper over forefoot and dorsum in Setaflex
- No seams at friction points
- Ultra lightweight anti-stumble flexible outsole



FLY

Style no. SA571

Size 35-41

Fitting 12 - Variable volume

Size chart B

- Flexpell® stretchable, breathable upper and patent leather with elastic strap
- No seams at friction points
- Flexible anti-stumble biomechanical outsole 20mm heel height
- Stabilising back counter with wide forepart for extra stability

ACTIVITY LINE

Designed to support walking and adapted physical activity. Advanced biomechanical technologies ensure support, control and safety at every step.

ADELLE SILVER VELCRO

Style no. SR717118

Size 35-41

Fitting 12.5

Size chart D

- Safeseam lining
- Soft Flexpell® upper over toes
- Biomechanical semi-rigid
- Extended back counter
- Super lightweight

FANCY BLACK

Style no. PA20171

Size 35-41

Fitting 12

Size chart C

- Safeseam lining
- Soft Flexpell® upper over toes
- Biomechanical semi-rigid
EVA Vibram® outsole
- Extended back counter
- Super lightweight





EVO MIDNIGHT VELCRO

Style no. PA20418

Size 35-41

Fitting 12.5

Size chart C

- Safeseam lining
- Super soft Flexpell® upper over toes
- Biomechanical semi-rigid Piuma® high abrasion resistance
- Super lightweight



POD X GREY

Style no. PA38774

Size 36-48

Fitting 11

Size chart J

- Knitted textile upper
- Safeseam construction
- Rearfoot external medial and lateral stabilisers for outstanding stability
- Super lightweight
- Semi-rigid high density D-lite sole
- Vegan

ACTIVITY LINE

Designed to support walking and adapted physical activity. Advanced biomechanical technologies ensure support, control and safety at every step.



SATURN BLUE

Style no. PA381380

Size 37-47

Fitting 11

Size chart J

- Knitted textile upper
- Safeseam construction
- Super lightweight
- Semi-rigid high density biomechanical sole
- Vegan



SATURN GREY

Style no. PA381371

Size 37-47

Fitting 11

Size chart J

- Knitted textile upper
- Safeseam construction
- Super lightweight
- Semi-rigid high density biomechanical sole
- Vegan



POD X BLACK

Style no. PA38771

Size 36-48

Fitting 11

Size chart J

- Knitted textile upper
- Safeseam construction
- Rearfoot external medial and lateral stabilisers for outstanding stability
- Super lightweight
- Semi-rigid high density D-lite sole
- Vegan



POD X NAVY

Style no. PA38770

Size 36-48

Fitting 11

Size chart J

- Knitted textile upper
- Safeseam construction
- Rearfoot external medial and lateral stabilisers for outstanding stability
- Super lightweight
- Semi-rigid high density D-lite sole
- Vegan

ACTIVITY LINE

Designed to support walking and adapted physical activity. Advanced biomechanical technologies ensure support, control and safety at every step.



FREEDOM EVO

Style no. PA34508

Size 35-46

Fitting 12

Size chart F

- Wide forepart for deformities, swelling, diabetes
- Water resistant nappa stretch
- Protected physical activity safeseam
- Semi-rigid biomechanical Timing Rocker® sole Piuma®
- Stabilising back counter

ITALA BLACK

Style no. PA4131

Size 35-41

Fitting 12

Size chart C

- Semi-rigid Timing Rocker® biomechanical outsole uses DCS (Dynamic Cross System®) technology
- Carbon fibre insert to curb excessive pronation/supination
- Non-collapsing outsole up to 140kg weight
- Soft Flexpell® safeseam upper



ACTIVITY ONE EVO VELCRO

Style no. PA30501

Size 39-47

Fitting 12

Size chart E

- Wide forepart for deformities, swelling or diabetes
- Protected physical activity safeseam
- Stabilising back counter



HERO EVO

Style no. PA30030

Size 39-46

Fitting 12

Size chart E

- Low stabilising back counter
- Biomechanical anti-stumbling Timing Rocker® Piuma® outsole
- Super lightweight

ACTIVITY LINE

Designed to support walking and adapted physical activity. Advanced biomechanical technologies ensure support, control and safety at every step.



SPORT DCS BLACK

Style no. PA37781

Size 39-47

Fitting 13

Size chart G

- Resistance against pronation and supination uses DCS (Dynamic Cross System®) technology
- Very high elastic response to help gait and reduce energy demand for heavier patients
- Stabilising back counter



SPORT DCS IRON

Style no. PA32780

Size 40-47

Fitting 14

Size chart L

- Textile/safeseam construction upper
- Safeseam construction
- Non collapsing outsole 140kg weight
- DCS rigid Timing Rocker® Vibram® sole



TURNER GREY

Style no. PA382478

Size 39-47

Fitting 12

Size chart H

- Perforated calf leather upper
- Safeseam construction
- D-lite DCS semi-rigid sole



TURNER GREY VELCRO

Style no. PA381479

Size 39-47

Fitting 12

Size chart H

- Perforated calf leather upper
- Safeseam construction
- D-lite DCS semi-rigid sole

ACTIVITY LINE



TURNER NAVY

Style no. PA381170

Size 39-47

Fitting 12

Size chart H

- Calf leather/textile upper
- Safeseam construction
- Leather qtrs/textile vamp
- D-lite DCS semi-rigid sole



WILD

Style no. PA365521

Size 39-47

Fitting 14.5

Size chart L

- Rigid biomechanical Timing Rocker® sole
- Extended stabilising heel counter

DIABETIC & RHEUMATIC

Footwear designed to support daily life with comfort, stability, and durability. Functional elegance with Podartis expertise.



GIOTTO

Style no. SD3001

Size 39-47

Fitting 11.5

Size chart H

- Semi-rigid biomechanical EVA sole with Vibram® outsole
- Heat mouldable upper
- Extended back counter
- Safeseam lining



WINDSOR BROWN SUEDE

Style no. SR30002

Size 39-46

Fitting 11.5

Size chart H

- Nubuck upper
- Safeseam construction
- Stabilising extended stiffener
- Flexible biomechanical

DIABETIC & RHEUMATIC

Footwear designed to support daily life with comfort, stability, and durability. Functional elegance with Podartis expertise.



Diabetic foot problem: severe deformities

X-DIAB 14

Style no. SD30091

Size 39-46

Fitting 14.5

Size chart G

- Biomechanical outsole and Vibram® energy saver technology
- Leather upper with Flexpell® stretchable vamp
- Stabilising back counter
- Removable insole allows insertion of custom-made insole
- Safeseam upper construction



Diabetic foot problem: severe deformities

X-DIAB 16

Style no. SD3471

Size 36-46

Fitting 16

Size chart K

- Biomechanical outsole and Vibram® energy saver technology
- Leather upper with Flexpell® stretchable vamp
- Stabilising back counter
- Removable insole allows insertion of custom-made insole
- Safeseam upper construction



BOTERO 11 VELCRO

Style no. SD370291

Size 39-46

Fitting 11.5

Size chart H

- Semi-rigid biomechanical EVA sole with Vibram® outsole
- Flexpell® stretchable upper to accommodate forefoot
- Extended back counter
- Safeseam lining



ENEZIA LAGUNA DARK GREY

Style no. SR680001

Size 36-41

Fitting 12

Size chart D

- Stretch textile/patent upper
- Safeseam construction
- Flexpell vamp
- EVA semi-rigid biomechanical rocker sole

ACUTE PHASE OFF-LOADING



TERAHEEL® HEEL PRESSURE OFF-LOADING


Style no. PO400

Size 35-46

Fitting Universal fit for left or right foot

Sold as singles

- Nylon upper
- Full opening with Velcro fastening
- Side stabilisers
- Customisable soft PU inlay
- Unique sole design for full protection and off-loading



STABIL D

Style no. PO700

Size 35-48

Fitting Universal fit for left or right foot

Sold as singles

- Off-loading and customisation
- Heat mouldable stabilising counter
- Removable insole to insert Modus or custom-made insole
- Rigid carbon fibre insole
- Rigid outsole
- Lateral stabilising removable whale bones



Stabil D + Modus insole
Lab Novel, 2007

Randomised double blind study on treatment in acute phase: Stabil-D vs TCC

Dr. Faglia, Dr. Caravaggi, Dr. Clerici, Dr. Sganzeroli, Dr. Curci, Dr. Simonetti,
Dr. Vailati, Dr. Sommalvico. Diabetic Foot Centres: IRCCS Multimedica -
Sesco S. Giovanni (Milan) Italy, Isturo Clinico Cita Studi (Milan) Italy.



TD-FOREFOOT & MIDFOOT

Style no. PO500

Size 35-48

Fitting Universal fit for left or right foot

Sold as singles

- No crutch: the biomechanical outsole enables the patient to walk safely and get back to activity sooner
- Preserves physiological posture: the outsole height is 22mm to minimise leg discrepancy
- Clinically tested for stability and effectiveness
- Removable insole to enable use with bespoke or our Modus insole for modular off-loading



MODUS INSOLE

Style No PL345

Size 35-48

Fitting Universal fit for left or right foot

Sold as singles

- Covered in dermocompatible DIAPOD, antibacterial and antifungal properties
- Ideal for plantar surface ulcer pressure relief

PODOWELL



ALADIN WIDE FIT WASHABLE STRETCH SHOE

Size 36-48

- A washable closed toe soft stretch vamp shoe, extra roomy with removable inlay
- Velcro fastening
- Polyurethane sole
- Black only



ALIX WIDE FIT WASHABLE SHOE

Size 36-46

- Microflex vamp, hand washable, wide fitting shoe with removable insole
- Velcro fastening
- Polyurethane sole
- Black only



ANDALOU WIDE FIT WASHABLE OPEN TOE

Size 36-46

- Open toe, soft textile upper with antibacterial finish and removable inlay for use with a captive Total Contact Insole
- Machine washable
- Wide fitting
- 6mm removable insole to replace with own
- Navy only



ATHOS WIDE FIT WASHABLE CLOSED TOE SHOE

Size 36-48

- A washable closed toe soft shoe, extra roomy with removable inlay
- Velcro fastening
- Polyurethane sole
- Black only

PODOWELL



ALTITUDE

Size 36-48

Fitting XL

- Ideal for use with oedematous or neuropathic feet, ulceration, incontinent users or swollen feet, and will accommodate bandaged feet
- Sanitized® antibacterial material limits development of bacteria and odours even after washing
- Anti-slip soles
- Machine washable at 40 degrees
- Fully opening for ease of donning
- Black only



XAVIER

Size 36-48

Fitting XXL

- Ideal for use with oedematous or neuropathic feet, ulceration, incontinent users or swollen feet, and will accommodate bandaged feet
- Microflex vamp, hand washable, wide fitting shoe with removable insole
- Velcro fastening
- Polyurethane sole
- Fully opening for ease of donning
- Grey only

VAROMED LYMPHOEDEMA SHOES



LYMPH VELCRO

Style no. 60330

Size 36-50 No size 49

- Velcro fastening, the vamp/tongue can be removed to produce a closed back sandal, elastic insert at the heel seat and a dual fastening system, 2 straps over the instep and a lace to secure around the leg
- Designed to accommodate bandaged legs and feet with a wide tongue, a lightweight sole and cushioned insole
- Black only



LYMPH LACE

Style no. 60310

Size 36-50 No size 49

- Lace fastening
- Designed to accommodate bandaged legs and feet with a wide tongue, a lightweight sole and cushioned insole
- Black only

PODARTIS FOOTWEAR SIZING GUIDE

CHART A												
SIZES	36	37	38	39	40	41	42	43	44	45	46	47
JOINT WIDTH	90	92	94	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	243	248	253	258	263	268	273	278	283	288	293	298
TOE BOX HEIGHT	32	33	34	35	36	37	38	39	40	41	42	43

PRODUCT RANGE
DEAMBULO L
DEAMBULO X

CHART B							
SIZES	35	36	37	38	39	40	41
JOINT WIDTH	87	89	90	92	93	95	96
JOINT CIRCUMFERENCE	238	241	246	253	258	261	268
TOE BOX HEIGHT	31	32	33	34	35	36	37

PRODUCT RANGE
FLY
DEAMBULO SPORT

CHART C							
SIZES	35	36	37	38	39	40	41
JOINT WIDTH	87	89	90	92	93	95	96
JOINT CIRCUMFERENCE	241	246	251	256	261	266	271
TOE BOX HEIGHT	36	37	38	39	40	41	42

PRODUCT RANGE
FANCY
ITALA
EVO MIDNIGHT

CHART D							
SIZES	35	36	37	38	39	40	41
JOINT WIDTH	85	87	88	90	91	93	94
JOINT CIRCUMFERENCE	241	246	251	256	261	266	271
TOE BOX HEIGHT	33	34	35	36	37	38	39

PRODUCT RANGE
ADELLE VELCRO
VENEZIA LAGUNA

CHART E									
SIZES	39	40	41	42	43	44	45	46	47
JOINT WIDTH	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	261	266	271	276	281	286	291	296	301
TOE BOX HEIGHT	34	34	34	35	36	37	38	39	40

PRODUCT RANGE
HERO EVO
ONE EVO MAN

CHART F												
SIZES	35	36	37	38	39	40	41	42	43	44	45	46
JOINT WIDTH	87	89	91	93	96	97	99	100	102	103	105	106
JOINT CIRCUMFERENCE	241	246	251	256	261	266	271	276	281	286	291	296
TOE BOX HEIGHT	34	34	34	34	34	34	34	35	36	37	38	39

PRODUCT RANGE
FREEDOM EVO

CHART G									
SIZES	39	40	41	42	43	44	45	46	47
JOINT WIDTH	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	261	266	271	276	281	286	291	295	301
TOE BOX HEIGHT	36	37	38	39	40	41	42	43	44

PRODUCT RANGE
DCS SPORT BLACK
X-DIAB 14

CHART H											
SIZES	37	38	39	40	41	42	43	44	45	46	47
JOINT WIDTH	94	95	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	246	251	256	261	266	271	276	281	286	291	296
TOE BOX HEIGHT	32	33	34	35	36	37	38	39	40	41	42

PRODUCT RANGE
GIOTTO
BOTERO 11 VELCRO
TURNER
WINDSOR

PODARTIS FOOTWEAR SIZING GUIDE

CHART J														PRODUCT RANGE	
SIZES	36	37	38	39	40	41	42	43	44	45	46	47	48	SATURN BLUE	
JOINT WIDTH	93	94	95	96	97	99	100	102	103	105	106	108	110	SATURN GREY	
JOINT CIRCUMFERENCE	247	251	256	261	266	271	276	281	286	291	296	301	306	POD X	
TOE BOX HEIGHT	34	34	35	35	36	36	37	38	38	39	40	40	41		

CHART K												PRODUCT RANGE	
SIZES	36	37	38	39	40	41	42	43	44	45	46	X-DIAB 16	
JOINT WIDTH	93	95	96	98	99	101	102	104	105	107	108		
JOINT CIRCUMFERENCE	262	267	272	277	282	287	292	297	302	307	312		
TOE BOX HEIGHT	39	40	41	42	43	44	45	46	47	48	49		

CHART L										PRODUCT RANGE	
SIZES	39	40	41	42	43	44	45	46	47	DCS SPORT IRON	
JOINT WIDTH	96	97	99	100	102	103	105	106	108	WILD	
JOINT CIRCUMFERENCE	261	266	271	276	281	286	291	296	301		
TOE BOX HEIGHT	35	36	36	37	38	38	39	40	40		

SIZE CHART FOR LYMPH SHOE/SANDAL ONLY

SIZE	36	37	38	39	40	41	42	43	44	45	46	47	48	50
FOOT LENGTH	239	246	253	259	266	273	279	286	293	299	306	313	319	333
JOINT WIDTH	91	93	95	97	99	101	103	105	107	109	111	113	115	119
HEEL WIDTH	60.5	62	63.5	65	66.5	68	69.5	71	72.5	74	75.5	77	78.5	81.5
TOE DEPTH	28.5	29	29.5	30	30.5	31	31.5	32	32.5	33	33.5	34	34.5	35.5
JOINT CIRC.	270	275	279	284	288	293	297	302	307	311	315	319	323	331
INSTEP CIRC.	280	285	289	294	298	303	307	312	317	321	325	329	333	341

reed

by eqwal.

SPRING 2026

Edition 5

www.reedmedical.co.uk

FOLLOW US ON



Scan me!

