

STEP INTO YOUR COMFORT ZONE

FOOTWEAR CATALOGUE



ISSUE 8





WELCOME TO THE 8TH EDITION OF THE REED MEDICAL THERAPEUTIC FOOTWEAR CATALOGUE.

WE NOW OFFER A FURTHER LAST SHAPE FOR THE AT RISK DIABETIC FOOT, THE **LOW TOE LAST** WHICH, THROUGH ITS INNOVATIVE DESIGN, ASSISTS IN THE REDUCTION OF PRESSURE UNDER THE METATARSAL HEADS.

Safe Seam® technology is now the proven standard in all of our diabetic footwear together with our unique “**Superplush Lining**” New designs in all categories of footwear are now available reflecting current high street trends.

We have increased choice in our popular leisure footwear range offering stylish new **TRAINERS** in a range of colours and new stretch materials exclusive to REED Medical.

Our well accepted **LYMPH** shoe is now available in Velcro as well as Lace fastening and can now be ordered in a modular specification.

Sales of our **PODARTIS** range continues to grow as the evidence based choice for ulcer management. The ever popular **PODOWELL** range has now been increased with a new style and more colour choice.

To complement our range of footwear we now offer a full range of Total Contact Insoles in Multi- Density EVA, Polypropylene , Polyurethane and Carbon Fibre, all with a choice of top covers and postings.

If you have any technical queries or would like any further information on any product in this catalogue please call us on 01254 503333.





LADIES - Classic Toe Shape

LADIES Classic



DIABETIC Collection

ALL STYLES IN THIS CATALOGUE
CAN BE ORDERED IN DIABETIC
SPECIFICATION



DIABETIC SPECIALIST - SHOWN IN BLACK

For those practitioners seeking less flexible footwear for the at risk diabetic foot we offer the SPECIALIST. Manufactured in "ARMY CALF" and with a moulded carbon fibre through stiffener, this footwear offers a more rigid environment in the treatment of the Diabetic foot. Superplush lining, through wedged soles and rocker with a sach heel completes the standard specification. Available in both Classic and Nature toe designs and in all widths and depths.

STANDARD SPECIALIST SPECIFICATION

- SAFESEAM® TECHNOLOGY
- SUPERPLUSH LINING
- RIM TOE PUFFS
- ARMY CALF
- CLASSIC, NATURE, LOW TOE ,
NOVA LAST SELECTION
- THROUGH WEDGE ROCKER SOLE
- CARBON FIBRE STIFFENERS
- SACH HEELS

OPTIONS

- MODULAR SPECIFICATION
- CHARCOT LAST ALLOWANCE
- TOTAL CONTACT INSOLES

BETH T-BAR
SHOWN IN BLACK LEATHER

A ladies Brogue T-Bar sandal
with single velcro fastening



BELLOWS
SHOWN IN TAUPE LEATHER

Extended opening bellows
tongue easily accommodates
variable oedema

DELICIA
SHOWN IN BURGUNDY
LEATHER /NUBUCK

Leather mudguard with
nubuck apron and quarters





FLEXTOP SHOE
SHOWN IN BLACK
LEATHER AND
NEOPRENE

Leather mudguard with
stretch neoprene vamp

RHIANNON
SHOWN IN MID
BROWN LEATHER

Seam free flower
punched vamp



ROSIE
SHOWN IN BROWN
LEATHER/NUBUCK

Extended opening style
leather with nubuck facing

KIERA
SHOWN IN
NAVY LEATHER

Normal opening wholecut
safe seam upper



THERAPY
SHOWN IN
BLACK LEATHER

Fully extended opening



CARINA
SHOWN IN
TAN NUBUCK

Normal opening shoe with
flower punch apron





FULLFLEX SHOE
SHOWN IN
TAUPE LEATHER AND
NEOPRENE

Stretch neoprene vamp
and soft padded collar

KERRY
SHOWN IN
ALABASTER LEATHER

Apron front with
basket weave



OPTIONS

BURGUNDY	
BROWN	
NAVY	
GREY	
TAUPE	
ALABASTER	

SALLY
SHOWN IN
BLACK LEATHER

Extended opening boot with
soft padded collar



PENNY
SHOWN IN
MID BROWN LEATHER

Flower punch apron with
mudguard and soft padded
collar



LOUISE
SHOWN IN
BLACK LEATHER

Leather boot with
seamfree vamp





ADIE
SHOWN IN
BROWN BRUISER

Rugged boot with plain front and cable stitch detail. Shown with optional commando sole.



CLAIRE
SHOWN IN
BLACK LEATHER

Leather plain vamp boot with nubuck trim



REBECCA
SHOWN IN
NAVY LEATHER

Seam free flower punched vamp

FLEXTOP BOOT
SHOWN IN
BROWN LEATHER AND
NEOPRENE

Neoprene apron with leather
mudguard



KIRSTEN
SHOWN IN
BROWN BRUISER

Rugged boot with wholecut
apron



REED CUSTOM INSOLES

Technology has always been an integral part of our business at Reed Medical, whether it is last carving or upper design we have always strived to create the perfect blend of both traditional and modern methods to create the high quality products.

It is with great excitement that we continue this tradition by announcing the CAD CAM insole range 'REED CUSTOM'. Available in Multi Density EVA, Polyurethane, Polypropylene and Carbon Fibre with multi choice coverings.

Once your orthotics are milled, they are expertly hand finished by a skilled technician to your exact requirements ensuring quality and perfect fit every time.



LADIES LEISURE - Classic Toe shape

LADIES Classic - LEISURE



DIABETIC HOUSE SHOE Collection

SEE PAGE 48-49
FOR FULL RANGE

Ideal for the 'At Risk'
patient who requires
a protective shoe in
the home



Continuing therapy in
the home mimics the
outdoor shoe using
indoor fabrics



HEATHER - SHOWN IN GREY

HEATHER -

FULLY CUSTOMISABLE

New water resistant airnet trainer, wholecut, with suede and leather laid on panels for a fully seam free construction and combined with our new absorber foam backed lining makes this product a comfortable, fashionable alternative to leather.



PINK



PURPLE



NAVY



BLACK



DENISE SHOWN IN WHITE SOFTEE/ BLUE PATENT FLASH

Soft leather plush lined
trainer with contrast flash





SERENA
SHOWN IN
WHITE SOFTIE/
PURPLE FLASH

Apron and mudguard with
contrast flash

DAISY
SHOWN IN
GREY

An extended opening
trainer with stretch suede
upper combination and full
stretch fitting. Ideal for the
sensitive foot.



VENUS
SHOWN IN
WHITE SOFTIE/
SILVER PATENT FLASH

Decorative stitch detail vamp
in soft leather





HELENA
SHOWN IN
BLACK/CERISE FLASH

Extended opening trainer
with contrast flash detail



PAULA
SHOWN IN
WHITE/ SOFT
BLUE FLASH

Extended opening trainer
boot with contrast flash



MARY
SHOWN IN
WHITE/ RED FLASH

Soft leather trainer boot with
contrast flash and back tab



HIDDEN DEPTHS - Nova Toe Shape

LADIES - Hidden Depths



SOLACE Collection

SEE PAGE 23-24
FOR FULL RANGE



THE SINGLE SOLUTION FOR THE RHEUMATOID FOOT

- Feminine styling
- Lightweight ultra soft leathers in pearl and satin finishes
 - Strong heel stiffeners for secure hold
 - Leather stacked heels or wedge units



VIOLET - WHITE LEATHER WITH
EMBROIDERED DETAIL

Flower embroidered detail shoe with PU sole

EMBROIDERY AVAILABLE ON ANY COLOUR

SUSAN
SHOWN IN
COPPER SATIN

Normal opening



ELLIE
SHOWN IN
NAVY SATIN

Extended opening shoe for
ease of donning



JANE
SHOWN IN
BURGUNDY LEATHER

Normal opening
Gibson shoe





NOVAFLEX
SHOWN IN
BLACK LEATHER/
NEOPRENE

Normal opening stretch
neoprene apron



BUNIONFLEX
SHOWN IN
BLACK LEATHER/
NEOPRENE

Stretch neoprene at the joint
to accommodate both
medial and lateral bunions



GILL
SHOWN IN
BURGUNDY LEATHER /
NUBUCK

Extended opening style with
nubuck detail facing

FLEUR
SHOWN IN
CHARCOAL SATIN

Delicate punch vamp & facing detail



ZOE
SHOWN IN
RED PEARL

New improved strap design on this popular style



VIOLET
SHOWN IN
BLACK LEATHER WITH
EMBROIDERED DETAIL

Flower embroidered detail





SARAH
SHOWN IN
BROWN NUBUCK/
GOLD

Casual bowling shoe with
contrast flash detail



JUDY
SHOWN IN
BLACK PEARL

Full toe T- bar sandal with
reverse pull velcro



LILY
SHOWN IN
BROWN NUBUCK

Normal opening shoe
with heavy decorative
stitch quarters

ERIN
SHOWN IN
BLACK LEATHER

Wholecut vamp with
collar detail



MEGAN
SHOWN IN
BLACK LEATHER

Plain vamp extended
opening boot



REED CUSTOM INSOLES

Technology has always been an integral part of our business at Reed Medical whether it is last carving or upper design we have always strived to create the perfect blend of both traditional and modern methods to create the high quality products.

It is with great excitement that we continue this tradition by announcing the CAD CAM insole range '**REED CUSTOM**'. Available in Multi Density EVA, Polyurethane, Polypropylene and Carbon Fibre with multi choice coverings.

Once your orthotics are milled, they are expertly hand finished by a skilled technician to your exact requirements ensuring quality and perfect fit every time.





LADIES SOLACE

LADIES Solace





SHARON
SHOWN IN
BEIGE LEATHER/GOLD &
SILVER SUEDE

An extended opening casual shoe with stud or diamante decoration with feminine locking grosgrain ribbon lace.

CHRISTINE
SHOWN IN NAVY SATIN

An extended opening casual shoe available in lace or our new more delicate 20mm velcro strap. This shoe can be made in any of our lovely new leather colours and finished in a light Eva wedge sole unit.



TRACEY
SHOWN IN WHITE
PEARL/WHITE SILVER SUEDE

Soft & comfortable leisure shoe made in white pearlescent leather with silver grosgrain ribbon lace.

LINDA
SHOWN IN
BLACK PEARL

A classic style brogue with extended opening for easy donning made in any of our soft pearlised leathers a leather stacked heel and flexible slip resistant rubber outsole and secure locking grosgrain ribbon lace.



CAROL
SHOWN IN
BEIGE PEARL

Attractive brogue punch gillie shoe with ribbon or lace fastening this style looks great and feminine with a dress or skirt. The leather stacked heel completes the look and looks equally as nice in all the new feminine look leathers.



JEAN
SHOWN IN BURGUNDY SATIN

An extended opening 6 eyelet shoe with a grosgrain ribbon lace and Eva wedge sole unit, choose the velcro option and the narrow straps enhance the delicate styling.





MENS - Classic Toe Shape

MENS Classic




OFFLOAD Collection

SEE PAGE 50-53
FOR FULL RANGE



**OFFLOADING FOR THE
AT RISK FOOT BY PODARTIS**

- Clinically tested
- Evidence based



GEORGE STRETCH BROGUE
SHOWN IN BROWN/ BROWN

NEW water resistant stretch fabric for ease of donning and swelling accommodation

THOR
SHOWN IN
BROWN BRUISER

Seamfree vamp with heavy
stitch detail on quarters



AUSTEN
SHOWN IN
BLACK LEATHER

Extended opening plain
vamp shoe



CALLUM
SHOWN IN
BLACK LEATHER

Normal opening plain
vamp shoe





THOMAS
SHOWN IN
MID BROWN LEATHER

Traditional brogue shoe




CONNOR
SHOWN IN
BROWN LEATHER

Normal opening heavy
stitch pintuck vamp



THERAPY
SHOWN IN
BLACK LEATHER

Fully extended opening



FLEXTOP SHOE
SHOWN IN
BLACK LEATHER /
NEOPRENE

Leather mudguard with
stretch neoprene vamp



BELLOWS
SHOWN IN
TAUPE LEATHER

Extended opening bellows
tongue easily accommodates
varying oedema



FULLFLEX SHOE
SHOWN IN
BLACK LEATHER /
NEOPRENE

Stretch neoprene vamp
and soft padded collar





ANDREW
SHOWN IN
BROWN LEATHER

Plain front wholecut
safeseam construction
with soft plush lining

REED CUSTOM INSOLES

Technology has always been an integral part of our business at Reed Medical whether it is last carving or upper design we have always strived to create the perfect blend of both traditional and modern methods to create the high quality products.

It is with great excitement that we continue this tradition by announcing the CAD CAM insole range '**REED CUSTOM**'. Available in Multi Density EVA, Polyurethane, Polypropylene and Carbon Fibre with multi choice coverings.

Once your orthotics are milled, they are expertly hand finished by a skilled technician to your exact requirements ensuring quality and perfect fit every time.



LORCAN
SHOWN IN
BLACK LEATHER

Normal opening plain
vamp boot



TYSON
MOCCA BRUISER

Rugged boot in bruiser
leather with contrast heavy
stitching on vamp & quarters

BEVAN
SHOWN IN
BLACK LEATHER

Classic brogue in a boot
styling





SAXON
SHOWN IN
BLACK LEATHER

Extended opening boot with
soft padded collar



FLEXTOP BOOT
SHOWN IN
BROWN LEATHER/
NEOPRENE

Neoprene apron with leather
mudguard



KYLE
SHOWN IN
BROWN BRUISER

Rugged boot with wholecut
apron



MENS LEISURE - Classic Toe Shape

MENS Classic - LEISURE



LYMPHOEDEMA Collection

SEE PAGE 54
FOR FULL RANGE



**AVAILABLE FROM STOCK
FOR FAST DELIVERY**

Can now also be made in modular or bespoke specification



MURRAY -
SHOWN IN BLUE



GREY



BLACK

MURRAY - SHOWN IN BLUE

New water resistant airnet trainer, wholecut, with suede and leather laid on panels for a fully seam free construction. Combined with our new absorber foam backed lining makes this product a comfortable, fashionable alternative to leather.

DALY
SHOWN IN
BLACK SOFTEE WITH
WHITE FLASH

Soft leather plush lined
trainer with contrast flash



COE
SHOWN IN
WHITE SOFTEE WITH
BLACK FLASH

Apron and mudguard with
contrast flash



BOLT
SHOWN IN
BLACK SOFTEE WITH
WHITE FLASH

Decorative stitch detail vamp
in soft leather





MURRAY
SHOWN IN
BLUE SUEDE/AIRNET

Extended opening trainer
boot with contrast flash

LENNOX
SHOWN IN
WHITE SOFTEE WITH
BLACK FLASH

Soft leather trainer boot



OVETT
SHOWN IN
WHITE SOFTEE WITH
BLACK FLASH

Extended opening trainer
boot with contrast flash





MENS HIDDEN DEPTH - Nova Toe Shape

MENS - HIDDEN DEPTH



HEAT MOULDABLE Collection


SEE PAGE 46-47
FOR FULL RANGE

Unique heat mouldable vamp for localised
customisation e.g. bunion width or clawed
toes etc.

Low temperature heat
gun or hair dryer

An insert is used
to bellow out the
material in specific
areas while heat is
applied

Once cooled the
material will hold its
shape until further heat
is applied





ATLANTA
SHOWN IN
BROWN NUBUCK

Extended opening style



TOLEDO
SHOWN IN
TAUPE NUBUCK

Full toe T-bar sandal with
reverse pull velcro



RENO
SHOWN IN
BLACK BRUISER

Seamfree vamp with
heavy stitch detail on
quarters

TUCSON
SHOWN IN
MID BROWN LEATHER

Gibson shoe



MONTANA
SHOWN IN
BROWN BRUISER

Casual shoe with contrast stitching



CARSON
SHOWN IN
MID BROWN LEATHER

Formal style with apron and mudguard detail





NOVAFLEX
BLACK LEATHER/
NEOPRENE

Normal opening stretch
neoprene apron



BUNIONFLEX
SHOWN IN
BLACK LEATHER/
NEOPRENE

Stretch neoprene at the joint
to accommodate both
medial and lateral bunions

VEGAS BROGUE
SHOWN IN
BLACK LEATHER

Classic brogue boot
with quarter punching &
safeseam construction



JOSEPH
SHOWN IN
BROWN NUBUCK

Plain front desert boot with
soft plush lining safeseam
construction



NEVADA
SHOWN IN
BLACK LEATHER

Normal opening boot
with wholecut vamp





MARK
SHOWN IN
BLACK LEATHER

Classic brogue desert boot with quarter punching & safeseam construction

REED CUSTOM INSOLES

Technology has always been an integral part of our business at Reed Medical, whether it is last carving or upper design we have always strived to create the perfect blend of both traditional and modern methods to create the high quality products.

It is with great excitement that we continue this tradition by announcing the CAD CAM insole range '**REED CUSTOM**'. Available in Multi Density EVA, Polyurethane, Polypropylene and Carbon Fibre with multi choice coverings.

Once your orthotics are milled, they are expertly hand finished by a skilled technician to your exact requirements ensuring quality and perfect fit every time.





REHABILITATION COLLECTION

- PODOWELL
- PODARTIS
- GAITOR INSOLES
- VANGUARD
- HEAT MOULDABLE
- LYMPHOEDEMA
- HOUSE SHOE

REHABILITATION Collection



GAITOR ARCH SUPPORTS

3/4 AND FULL LENGTH GAITOR ARCH SUPPORTS FOR SHOCK ABSORPTION, ARCH SUPPORT & COMFORT

Gaitor 3/4 and Full Length Arch Supports are designed to provide superior shock absorption, arch support and comfort to your feet. Incorporating a successfully patented waffle design which effectively duplicates nature's own shock absorbing system in your feet, the combination of modern materials provide a lightweight, yet highly effective shock absorber.

Revolutionary new technologies offer an ultra thin carbon fibre arch support which will provide long-lasting relief and comfort. This technologically advanced carbon material may be adjusted when necessary, offering a degree of customisation.

Instant relief for

- Painful fallen arches
- Painful heel conditions and heel spurs
- Sore, tired, burning, aching feet
- Shin splints
- Painful legs, ankles and knees
- Painful neck and back
- Arthritic joints

In addition Gaitors contain a combination of materials which keep your feet cooler and drier. Bacterial and fungal growth inhibitors have been added to increase sanitation and reduce unwanted odours.

Size Range

S 3-4

M 5-7

L 8-10

XL 11-14

Features

- Mouldable Support
- Cushioned Heel
- Plantar Fasciitis and Heel Spur Relief
- Stay Fresh Fabric





VANGUARD SHOE - BLACK

Fully opening shoe ideal over bandages and swelling



VANGUARD BOOT
BLACK

Fully opening boot ideal over bandages and swelling



NEOPRENE BOOT
BLACK/ BLACK LEATHER

Full stretch neoprene boot with heel stiffener for stability. Available in lace or velcro. Can be made to specific ankle and leg measurements

YEW
Captive Velcro
fastening



ASH
Available in
lace or Velcro
fastening



FABRIC SWELLIES - BLACK ONLY

- Heat mouldable vamp for specific bunion width and claw toe depth
- Suitable for both indoor and outdoor use
- Hand washable at 30°C
- Superplush anti-bacterial lining
- Suitable for orthotic adaptations
- Available in stock sizes 1-12 in narrow, standard, wide, deep and wide & deep
- Modular options
- Rigid rocker sole option
- Lightweight

ASPEN

Captive Velcro
fastening



ELM

Available in lace or
Velcro fastening



LEATHER LOOK SWELLIES - BLACK ONLY

- A full leather look option is now available with heat mouldable vamp
- Suitable for both indoor and outdoor use
- Hand washable at 30°C
- Superplush anti-bacterial lining
- Suitable for orthotic adaptations
- Available in stock sizes 1-12 in narrow, standard, wide, deep and wide & deep
- Modular options
- Rigid rocker sole option
- Lightweight

PISCES
BURGUNDY

Soft velveteen uppers
velour lining with a seam
free vamp



LIBRA
**AVAILABLE IN BEIGE
AND PINK**

Soft velveteen uppers and
fleece lining with a seam free
vamp



AQUARIUS
**AVAILABLE IN BEIGE
AND PINK**

Soft velveteen uppers
and fleece lining with a
seam free vamp






GEMINI BURGUNDY

Burgundy soft velveteen upper with soft velour lining

LEO AVAILABLE IN NAVY OR BURGUNDY CHECK

Classic check tweed with warm and soft velour lining



TAURUS AVAILABLE IN BURGUNDY OR NAVY CHECK

Open toe, soft textile upper with anti-bacterial finish and removable inlay for use with a captive Total Contact Insole



FOREFOOT OFFLOADER - PO100 SIZES 35-48

Offloading of ulcers, pressures
on the forefoot and ulcers
between toes



INDICATIONS

Ideal angle of ankle
dorsiflexion to stabilize the
walk and reduce the “step
effect”

Fully opening with triple
velcro closure

Removable EVA insole

8° degrees outsole with
ankle dorsiflexion to transfer
the pressure to the rear foot

Outsole height only 34mm

EVIDENCE BASED OFFLOADING

kPa

300.00

220.00

150.00


100.00

60.00

40.00

15.00

15.00





HEEL OFFLOADER PO400

SIZES S: 35-37, M: 38-40,
L: 41-43, XL: 44-46

Offloading of Heel



INDICATIONS

Complete offloading of the
rearfoot

Rear removable protection

Customizable offloading
insole in soft PU

Fully opening with 3 velcro
straps

Outsole tested to offload the
heel completely

EVIDENCE BASED OFFLOADING

kPa

300.00

220.00

150.00


100.00

60.00

40.00

15.00

15.00



FOREFOOT AND MIDFOOT OFFLOADER PO500

SIZES 35-48

Offloading of digital ulcers,
midfoot ulcers, lesions and
ulcers on the forefoot



INDICATIONS

No need to use crutches: the biomechanical outsole enables the patient to walk safely and to get back sooner to a normal life

Stabilising back counter

22mm outsole height

Fully opening 3 velcro strap fastening

Biomechanical rigid outsole

Preserves the physiological posture: The outsole height is 22mm to minimise the leg length discrepancy

Clinically tested for stability and effectiveness

EVIDENCE BASED OFFLOADING

kPa

300.00

220.00

150.00


100.00

60.00

40.00

15.00

15.00





STABILITY WALKER PO700 SIZES 35-46

Treatment of neuropathic and ischemic ulcers in the midfoot / forefoot
Charcot foot. Particularly recommended in the treatment of diabetic vascular foot

INDICATIONS

Offloading and customisation
Rear foot heat mouldable stabilising bones
Space for Custom / Modus insole
Rigid rocker outsole to aid stability on the midfoot stance
Tongue to protect the tibia and throat from friction



DUAL DENSITY 8 DEGREE NEGATIVE PITCH INSOLE - PO060

Transfers pressure from the forefoot to the rearfoot

SIZES 35/36 37/38 39/40 41/42 43/44 45/46 47/48

Universal fit left or right



MODUS INSOLE - PL345

Triple density insole with removable inserts for specific pressure reduction

Anti bacterial and anti fungal properties ideal for planter, surface, ulcer pressure relief

SIZES 35/36 37/38 39/40 41/42 43/44 45/46 47/48

Universal fit left or right

See accomodation
charts for sizes 36-50

**60330 LYMPH
SHOE VELCRO**
BLACK STRETCH VELOUR

SIZES 36-50

Upper micro/stretch
combined, elastic insert
at the heel seat, velcro
fastening



**60310 LYMPH
SHOE LACE**
BLACK STRETCH VELOUR

SIZES 36-50

Upper micro/stretch
combined, elastic insert
at the heel seat, lacing,
removable 3mm soft
foam insole

**60320 LYMPH
SANDAL**
BLACK STRETCH VELOUR

SIZES 36-50

Stretch Upper, velcro
fastening, soft foam insole.
Available with toe protector
(Sold seperatly 60321)





REHABILITATION PodoWell

REHABILITATION PodoWell



ALADIN
BLACK ONLY
SIZES 36-48

A washable shoe with stretch microflex vamp extra roomy ideal for people with Hallux Valgus (Bunion).



ATHOS
AVAILABLE IN NAVY
AND BLACK
SIZES 36-48

A washable closed toe soft shoe extra roomy with removable inlay



AJACCIO
BLACK ONLY
SIZES 36-46

Microflex vamp, hand washable, wide fitting shoe with removable insole





AWELL

AVAILABLE IN
BROWN & BEIGE
SIZES 36-46

A washable shoe with stretch microflex panels extra roomy ideal for people with Hallux Valgus (Bunion).



ANDALOU

NAVY ONLY
SIZES 36-46

Open toe, soft textile upper with anti-bacterial finish and removable inlay for use with a captive Total Contact Insole



MEASURING FOR REED FOOTWEAR



SHOE SIZE

Measure the length of the foot using your Reed measuring stick. Place the fixed jaw of the stick at the heel and measure to the longest toe. Hold the sliding jaw in place and read the corresponding size on the stick. Add 1 1/2 sizes if patient standing or 2 if sitting. This is the shoe size.

JOINT WIDTH

Open the jaw on your stick place either side of the joint, hold the jaw in place and read the measure.

TOE DEPTH

Place fixed jaw under toes and measure depth to highest toe, hold the jaw in place and read the measure.

JOINT CIRCUMFERENCE

Place your tape measure under the met heads and measure around the widest part of the foot.

INSTEP CIRCUMFERENCE

Place your tape measure under the waist, meeting at the navicular.

LEG CIRCUMFERENCE A: MEASURE (BOOTS ONLY)

Place your tape measure at 125mm height and measure around the leg

ANKLE CIRCUMFERENCE B: MEASURE (BOOTS ONLY)

Place your tape measure around the ankle bone to determine **B** measure

INSOLE ALLOWANCE

Consider the depth of insole to be used. All footwear has 9mm allowance excluding Lowline 3mm and Solace 6mm.

FITTING CHARTS


Refer to fitting chart to find the closest match for your measures. We can't reduce a last. If you require increased measures specify your modular adaptations.

All our footwear is supplied with 3x3mm inlays excluding our Lowline last which has 1x3mm inlay and Solace which has 1x6mm inlay.

The footwear is supplied with sole as catalogue if you require an alternative please specify when you order.

Standard leg measures on boots will be supplied unless custom measurements are specified. When supplying these please specify the height the measure was taken.

A: Circumference of leg at 125mm"
B: Circumference at ankle bone



LADIES ACCOMODATION CHART

LADIES NARROW LOWLINE

3mm Insole allowance 1x3mm footbed

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
1	77mm	206mm - 220mm	190mm - 200mm	28mm	54mm
2	79mm	213mm - 227mm	196mm - 206mm	29mm	55mm
3	81mm	219mm - 234mm	200mm - 212mm	30mm	56mm
4	83mm	227mm - 239mm	206mm - 220mm	31mm	57mm
5	85mm	232mm - 246mm	212mm - 227mm	32mm	58mm
6	87mm	235mm - 248mm	220mm - 232mm	33mm	59mm
7	89mm	241mm - 254mm	227mm - 239mm	34mm	60mm
8	91mm	248mm - 260mm	232mm - 246mm	35mm	62mm
9	93mm	255mm - 267mm	241mm - 253mm	36mm	63mm

LADIES STANDARD

9mm Insole allowance 3x3mm footbed

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
1	87mm	222mm - 231mm	209mm - 222mm	31mm	55mm
2	89mm	228mm - 235mm	216mm - 229mm	32mm	56mm
3	91mm	235mm - 248mm	222mm - 235mm	33mm	57mm
4	93mm	241mm - 254mm	228mm - 241mm	34mm	58mm
5	95mm	248mm - 260mm	235mm - 248mm	35mm	59mm
6	97mm	254mm - 266mm	241mm - 254mm	36mm	60mm
7	99mm	260mm - 273mm	248mm - 260mm	37mm	61mm
8	101mm	266mm - 279mm	254mm - 266mm	38mm	62mm
9	103mm	272mm - 286mm	260mm - 272mm	39mm	63mm

LADIES WIDE LOWLINE

3mm Insole allowance 1x3mm footbed

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
1	93mm	235mm - 248mm	216mm - 228mm	31mm	56mm
2	95mm	241mm - 254mm	222mm - 235mm	32mm	57mm
3	97mm	248mm - 260mm	228mm - 241mm	33mm	58mm
4	99mm	254mm - 266mm	235mm - 248mm	34mm	59mm
5	100mm	260mm - 273mm	241mm - 254mm	35mm	60mm
6	103mm	266mm - 279mm	248mm - 260mm	36mm	61mm
7	107mm	273mm - 287mm	254mm - 266mm	37mm	62mm
8	109mm	279mm - 293mm	260mm - 273mm	38mm	63mm
9	111mm	285mm - 299mm	266mm - 279mm	39mm	64mm

LADIES NARROW
9mm Insole allowance 3x3mm footbed

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
1	77mm	206mm - 220mm	190mm - 200mm	31mm	54mm
2	79mm	213mm - 227mm	196mm - 206mm	32mm	55mm
3	81mm	219mm - 234mm	200mm - 212mm	33mm	56mm
4	83mm	227mm - 239mm	206mm - 220mm	34mm	57mm
5	85mm	232mm - 246mm	212mm - 227mm	35mm	58mm
6	87mm	235mm - 248mm	220mm - 232mm	36mm	59mm
7	89mm	241mm - 254mm	227mm - 239mm	37mm	60mm
8	91mm	248mm - 260mm	232mm - 246mm	38mm	61mm
9	93mm	255mm - 268mm	237mm - 252mm	39mm	62mm

LADIES STANDARD LOWLINE
9mm Insole allowance 3x3mm footbed

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
1	87mm	222mm - 231mm	209mm - 222mm	28mm	55mm
2	89mm	228mm - 235mm	216mm - 229mm	29mm	56mm
3	91mm	235mm - 248mm	222mm - 235mm	30mm	57mm
4	93mm	241mm - 254mm	228mm - 241mm	31mm	58mm
5	95mm	248mm - 260mm	235mm - 248mm	32mm	59mm
6	97mm	254mm - 266mm	241mm - 254mm	33mm	60mm
7	99mm	260mm - 273mm	248mm - 260mm	34mm	61mm
8	101mm	266mm - 279mm	254mm - 266mm	35mm	62mm
9	103mm	272mm - 285mm	260mm - 272mm	36mm	63mm

LADIES WIDE
9mm Insole allowance 3x3mm footbed

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
1	93mm	235mm - 248mm	216mm - 228mm	31mm	56mm
2	95mm	241mm - 254mm	222mm - 235mm	32mm	57mm
3	97mm	248mm - 260mm	228mm - 241mm	33mm	58mm
4	99mm	254mm - 266mm	235mm - 248mm	34mm	59mm
5	101mm	260mm - 273mm	241mm - 254mm	35mm	60mm
6	103mm	266mm - 279mm	248mm - 260mm	36mm	61mm
7	105mm	273mm - 287mm	254mm - 266mm	37mm	62mm
8	107mm	279mm - 293mm	260mm - 273mm	38mm	63mm
9	109mm	285mm - 300mm	267mm - 279mm	39mm	64mm

LADIES DEEP 9mm Insole allowance 3x3mm footbed					
SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
1	87mm	235mm - 248mm	216mm - 228mm	40mm	55mm
2	89mm	241mm - 254mm	222mm - 235mm	41mm	56mm
3	91mm	248mm - 260mm	228mm - 241mm	42mm	57mm
4	93mm	254mm - 266mm	235mm - 248mm	43mm	58mm
5	95mm	260mm - 273mm	241mm - 254mm	44mm	59mm
6	97mm	266mm - 279mm	248mm - 260mm	45mm	60mm
7	99mm	273mm - 287mm	254mm - 266mm	46mm	61mm
8	101mm	279mm - 293mm	260mm - 273mm	47mm	62mm
9	103mm	285mm - 299mm	266mm - 279mm	48mm	63mm

LADIES WIDE & DEEP 9mm Insole allowance 3x3mm footbed					
SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
1	93mm	235mm - 248mm	216mm - 228mm	40mm	56mm
2	95mm	241mm - 254mm	222mm - 235mm	41mm	57mm
3	97mm	248mm - 260mm	228mm - 241mm	42mm	58mm
4	99mm	254mm - 266mm	235mm - 248mm	43mm	59mm
5	100mm	260mm - 273mm	241mm - 254mm	44mm	60mm
6	103mm	266mm - 279mm	248mm - 260mm	45mm	61mm
7	107mm	273mm - 287mm	254mm - 266mm	46mm	62mm
8	109mm	279mm - 293mm	260mm - 273mm	47mm	63mm
9	112mm	285mm - 299mm	266mm - 279mm	48mm	64mm

SIZE CHART FOR LYMPH SHOE/SANDAL ONLY														
SIZE	36	37	38	39	40	41	42	43	44	45	46	47	48	50
Foot Lenght	239	246	253	259	266	273	279	286	293	299	306	313	319	333
Joint Width	91	93	95	97	99	101	103	105	107	109	111	113	115	119
Heel Width	60.5	62	63.5	65	66.5	68	69.5	71	72.5	74	75.5	77	78.5	81.5
Toe Depth	28.5	29	29.5	30	30.5	31	31.5	32	32.5	33	33.5	34	34,5	35.5
Joint Circ.	270	275	279	284	288	293	297	302	307	311	315	319	323	331
Instep Circ.	280	285	289	294	298	303	307	312	317	321	325	329	333	341

MENS ACCOMODATION CHART

MENS NARROW LOWLINE

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
5	85mm	224mm - 234mm	222mm - 232mm	32mm	58mm
6	88mm	230mm - 240mm	228mm - 238mm	33mm	59mm
7	89mm	236mm - 246mm	234mm - 244mm	34mm	60mm
8	91mm	242mm - 252mm	240mm - 250mm	35mm	61mm
9	93mm	249mm - 259mm	247mm - 257mm	36mm	62mm
10	95mm	255mm - 265mm	253mm - 263mm	37mm	63mm
11	97mm	263mm - 273mm	261mm - 271mm	38mm	64mm
12	99mm	269mm - 279mm	267mm - 277mm	39mm	65mm

MENS STANDARD LOWLINE

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
5	94mm	236mm - 246mm	234mm - 244mm	32mm	61mm
6	95mm	242mm - 252mm	240mm - 250mm	33mm	62mm
7	97mm	249mm - 259mm	246mm - 256mm	34mm	63mm
8	99mm	255mm - 265mm	253mm - 263mm	35mm	64mm
9	102mm	263mm - 273mm	261mm - 271mm	36mm	65mm
10	105mm	269mm - 279mm	267mm - 277mm	37mm	66mm
11	108mm	275mm - 285mm	273mm - 283mm	38mm	67mm
12	112mm	282mm - 292mm	280mm - 290mm	39mm	68mm

MENS WIDE LOWLINE

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
5	99mm	246mm - 256mm	238mm - 248mm	32mm	61mm
6	102mm	253mm - 263mm	245mm - 255mm	33mm	62mm
7	105mm	259mm - 269mm	252mm - 262mm	34mm	63mm
8	108mm	265mm - 275mm	259mm - 269mm	35mm	64mm
9	112mm	272mm - 282mm	266mm - 276mm	36mm	65mm
10	115mm	278mm - 288mm	272mm - 282mm	37mm	66mm
11	118mm	286mm - 296mm	278mm - 288mm	38mm	67mm
12	120mm	294mm - 302mm	284mm - 294mm	39mm	68mm



MENS NARROW

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
5	85mm	235mm - 248mm	228mm - 241mm	35mm	58mm
6	88mm	241mm - 254mm	235mm - 248mm	36mm	59mm
7	89mm	248mm - 260mm	241mm - 254mm	37mm	60mm
8	91mm	254mm - 266mm	248mm - 260mm	38mm	61mm
9	93mm	260mm - 273mm	254mm - 266mm	39mm	62mm
10	95mm	266mm - 279mm	260mm - 273mm	40mm	63mm
11	97mm	273mm - 287mm	266mm - 279mm	41mm	64mm
12	99mm	279mm - 293mm	273mm - 287mm	42mm	65mm

MENS STANDARD

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
5	94mm	248mm - 260mm	241mm - 254mm	35mm	61mm
6	95mm	254mm - 266mm	248mm - 260mm	36mm	62mm
7	97mm	260mm - 273mm	254mm - 266mm	37mm	63mm
8	99mm	266mm - 279mm	260mm - 273mm	38mm	64mm
9	102mm	273mm - 287mm	266mm - 279mm	39mm	65mm
10	105mm	279mm - 293mm	273mm - 287mm	40mm	66mm
11	108mm	287mm - 299mm	279mm - 293mm	41mm	67mm
12	112mm	293mm - 306mm	287mm - 299mm	42mm	68mm

MENS WIDE

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
5	99mm	258mm - 270mm	248mm - 254mm	35mm	61mm
6	102mm	264mm - 277mm	252mm - 260mm	36mm	62mm
7	105mm	270mm - 283mm	254mm - 264mm	37mm	63mm
8	108mm	277mm - 290mm	264mm - 270mm	38mm	64mm
9	112mm	283mm - 296mm	270mm - 277mm	39mm	65mm
10	115mm	290mm - 302mm	277mm - 283mm	40mm	66mm
11	118mm	299mm - 311mm	285mm - 295mm	41mm	67mm
12	120mm	308mm - 320mm	295mm - 301mm	42mm	68mm

MENS DEEP

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
5	94mm	266mm - 279mm	260mm - 273mm	44mm	61mm
6	95mm	273mm - 287mm	266mm - 279mm	45mm	62mm
7	97mm	279mm - 293mm	273mm - 287mm	46mm	63mm
8	99mm	287mm - 299mm	279mm - 293mm	47mm	64mm
9	102mm	293mm - 306mm	287mm - 299mm	48mm	65mm
10	105mm	299mm - 312mm	293mm - 306mm	49mm	66mm
11	108mm	306mm - 318mm	299mm - 312mm	50mm	67mm
12	112mm	308mm - 320mm	306mm - 314mm	51mm	68mm

MENS WIDE & DEEP

SIZE	JOINT WIDTH	JOINT CIRCUMFERENCE	INSTEP CIRCUMFERENCE	TOE BOX DEPTHS	HEEL WIDTH
5	99mm	266mm - 280mm	259mm - 273mm	44mm	61mm
6	102mm	274mm - 288mm	265mm - 279mm	45mm	62mm
7	105mm	282mm - 296mm	272mm - 286mm	46mm	63mm
8	108mm	290mm - 304mm	279mm - 293mm	47mm	64mm
9	112mm	298mm - 312mm	285mm - 299mm	48mm	65mm
10	115mm	306mm - 320mm	292mm - 306mm	49mm	66mm
11	118mm	314mm - 328mm	298mm - 312mm	50mm	67mm
12	120mm	322mm - 336mm	304mm - 318mm	51mm	68mm


SIZE CHART FOR LYMPH
SHOE/SANDAL ONLY

SIZE	36	37	38	39	40	41	42	43	44	45	46	47	48	50
Foot Lenght	239	246	253	259	266	273	279	286	293	299	306	313	319	333
Joint Width	91	93	95	97	99	101	103	105	107	109	111	113	115	119
Heel Width	60.5	62	63.5	65	66.5	68	69.5	71	72.5	74	75.5	77	78.5	81.5
Toe Depth	28.5	29	29.5	30	30.5	31	31.5	32	32.5	33	33.5	34	34.5	35.5
Joint Circ.	270	275	279	284	288	293	297	302	307	311	315	319	323	331
Instep Circ.	280	285	289	294	298	303	307	312	317	321	325	329	333	341

LAST ADAPTATION INFORMATION


MODULAR LAST OPTIONS

Additions are available in multiples of 3mm any adaption over 6mm will require a Bespoke pattern.



LAST SHAPES

N.B. Narrow heel is available on our Classic and Nova last only, ask for Narrow heels when ordering your footwear



This last design has a reduction of 6mm through depth from the regular last to accommodate those patients whose prescription demands a narrow, standard or wide fitting, yet do not require the extra depth associated with the regular last shape.

We now offer a further last shape for the at risk diabetic foot, the LOW TOE LAST which, through its innovative design, assists in the reduction of pressure under the metatarsal heads.

LEATHER

	BLACK
	GREY
	TAUPE
	BROWN
	BROWN ANTIQUE BRUISER
	MOCKA BRUISER
	TAN
	NAVY
	RED
	BURGUNDY
	ALABASTER
	WHITE

PATENT

	BLACK
	SILVER
	NAVY BLUE
	RED
	CERISE
	PINK
	LILAC
	AUBERGINE
	WHITE

NUBUCK

	BLACK
	GREY
	TAUPE
	BROWN
	CHOCOLATE
	TAN
	NAVY
	SOFT BLUE
	BURGUNDY
	PINK
	LILAC
	FUSCIA
	TERRACOTTA
	HONEY
	SAGE
	PISTACHIO

SUEDE

	GREY
	BROWN
	NAVY

SATIN

	BURGUNDY
	NAVY
	CHARCOAL
	COPPER
	VIOLET

PEARL

	BLACK
	WHITE
	RED
	MAGENTA
	BEIGE
	SOFT BLUE
	SOFT PINK

SOLING OPTIONS



PYRAMID EVA



VIBRAM COMMANDO



COMMANDO BLACK



HIDDEN DEPTH PU BLACK



COMMANDO NATURAL



TCR OIL RESISTANT



TWS COMMANDO



DOG TOOTH RUBBER MIX

SOLE PATTERN OPTIONS



ULLI EVA



WELLAN



RIBBED RUBBER



HARD RUBBER



LIGHTWEIGHT RIBBED



CARACOLAS RUBBER



TYRE TRACK RUBBER

TERMS AND CONDITIONS

All of our footwear can be made with a diabetic specification of softlee leather, rim toe puff & plush lining.

UK despatch is next day service. Next day timed delivery is available but will incur extra cost

All adaptations are carried out above the feather edge so as not to alter the bottom pattern so that a standard unit will fit. Should bottom pattern alterations be necessary this will then become bespoke and may incur further costs.

All footwear is supplied with 3 x 3mm inlays. Excluding Lowline which is supplied with 1 x 3mm and Solace 1 x 6mm

Standard sole & heels and through wedge soles are made from EVA.

Colour reproduction - We make every effort to ensure the colour you see in our catalogue and on our website is accurate, however, colours can be different in photographs and shades can vary from batch to batch and over the life of our catalogue.

Only Callum, Austin, Keira and Therapy made in Black, Brown or Navy which are made on our Reeco classic last is returnable without any alteration / adaption.

1st trial fitting is free of charge and subsequent fittings may be charged for. Accurate measurements are essential when placing your order as re-lasting the upper and new uppers etc following trial fitting will be charged for. An invoice will be raised for all orders not returned for completion with 90 days.

Split sizes, split fittings, alterations to lasts and special orders will be charged in accordance with the current price list and are non-returnable.

We only use return Velcro fastening straps on all our shoes to enable greater adjustability. Lay on velcro fastening is available as an option but is then non- returnable.

Returnable footwear must be in a resaleable condition and received by us within 1 calendar month following despatch for a full refund.

All shelf stock footwear - Podowell, Podartis and Lymph are returnable with a full refund unworn and in the original packaging.

Documentation relating to the confirmation of orders should be clearly marked 'confirmation' to avoid duplication. Duplicated orders through originators error will be chargeable.

The company reserves the right to alter specification or withdraw styles without prior notice. Specification changes are only undertaken when product improvements are identified or a component material is no longer available.

All goods remain the property of Reed Medical Ltd until paid for in full.

RETURNS

Both men's & ladies boots are non-returnable.

Single footwear is non-returnable.

Worn footwear is non-returnable.

Footwear made on our Lowline , Nova and Nature last is non-returnable.

Footwear ordered in a non- standard colour or combination is non-returnable.

WEBSITE

All our catalogues and brochures are available for download from our website www.reedmedical.co.uk our website is a live document and we constantly feed it with information regarding all of our new styles, current news items etc. Please visit it regularly or follow us on facebook and you will receive updates/ information regularly and as it happens.



REED MEDICAL LTD.

Thwaites Close
Shadsworth Business Park
Blackburn, Lancashire
United Kingdom BB1 2QQ

Telephone: 01254 503333
Fax: 01254 503360
E-mail: admin@reedmedical.co.uk

www.reedmedical.co.uk

MENS ACCOMODATION CHART

LADIES ACCOMODATION CHART