




REED Medical Footwear Catalogue **2018**



Kids Collection

Our shoes have been individually hand crafted with great care.



Made in Britain



Anything you
want, just ask!



REED



Kids Collection



For over 60 years Reed has become known in the industry as a manufacturer of Quality footwear for Children and Adults.

Based in Lancashire our highly experienced workforce make stylish quality footwear using the best leathers, suedes, patents, nubucks and fabrics. Using high street styling as a design tool we create a range that doesn't stand out in the crowd!

We have listened to our customers and this year we launch our Kidztuff® range a more robust construction with:

- Medially extended heavy duty heel stiffener
- Robust upper leathers
- Full leather lining
- Rigid rubber sole unit
- Increased backheight for stability

Choose from our Fashionable:

- Tiny toes
- Kids collection
- School collection
- Kidztuff® collection
- Stability & Toewalkers
- Straight/Reverse and AFO
- Ponseti boots & Dobbs bars
- New INSOLE selection

Require further durability? Add panels, mudguard or outside counter in our new Anti Scuff material which is sewn on to ensure our styles keep their cosmesis.

We have a dedicated customer service and technical team ready to answer any questions you might have.

Don't see what you want?
Ring our team and "just ask".
We are the manufacturer and can meet all your requirements.



Contents

Tiny Toes 15-20	4
Kids Collection 21 - 40	6
Insoles	17
School Collection 21 - 40	18
Kidztuff Collection	22
AFO Collection	26
Straight & Reverse	28
Stability & Toe Walkers	29
Ponsetti Boots	30
Dobbs Bars	31
Colour Chart	32
Embroidery Options	33
Accommodation Chart	34





Tiny Toes 15-20



Baby Vyper

Extended opening trainer boot with contrast flash.
 Shown in white leather & soft pink with plush lining.



Bubble

Open to toe boot with punch detail with animal print collar & outside back.
 Shown in pistachio nubuck.



Raven

Extended opening plain boot & plush lining.
 Shown in white & soft pink.



Squeak/Teddy

Open to toe boot with teddy bear tongue & plush lining.
 Shown in navy nubuck & beige leather.



Squeak/Unicorn

Extended opening boot with unicorn embroidery.
 Shown in pink leather.



Choose your own design images on page 25

Whatever you want... just ask



Kids Collection 21 - 40

Ruby

Court shoe style sandal with heart detail. Velcro fasten.


Shown in black leather.



Opal

Court shoe with single rose embroidery.

Shown in black patent.



Embroidery Options

Sapphire

Court shoe style sandal. Plain front. Velcro fasten.

Shown in red patent.



Holly

Normal opening sandal with flower detail. Double strap velcro fasten.

Shown in white patent and black leather.



Jasmine

Normal opening sandal. Double strap velcro fasten with flower detail.

Shown in violet patent.



Sophie

Normal opening double return velcro fastening.

Shown in cobalt patent & navy nubuck heart detail.





Alex

3 Strap velcro sandal with padded collar.

Shown in navy leather & leather lining.



Moon

2 Strap velcro sandal with padded collar.

Shown in brown bruiser.



Dodger

Normal opening boot.

Shown in grey patent and black leather.



Raven

Extended opening boot.

Shown in brown bruiser and black leather.



Dune

3 Strap velcro sandal with padded collar.

Shown in pistachio/cream nubuck.



Neurotec modular foot
orthotic system - see page 17



Whatever you want... just ask 9



Blossom

Normal opening with daisy quarters, please see colour chart below. Commando sole as standard.

Shown in red leather, honey nubuck and black leather.



Available in black or red daisy. Red, Black, Cerise, pink, Pistachio, Beige Animal print.



Josie

Normal opening desert boot style with multi coloured thread detail. Commando sole as standard.

Shown in purple nubuck.



Cougar

Extended opening rugged boot with leather contrast collar. Commando sole as standard.

Shown in honey nubuck.



Millie

Normal opening boot with check nubuck or animal print collar.

Shown in cerise leather and red patent.



Star

Extended opening boot with animal print quarters and ski hook fastening. Commando sole as standard.

Shown in soft pink leather & animal print.



Whatever you want... just ask ||

Retro

Popular trainer shoe.

Shown in red nubuck, navy leather & natural leather all with contrast flash.

Choose your colour combination...



Spyder

Trainer boot with padded collar.


Shown in black leather with white flash and white leather with silver flash.



Adder

Water resistant Airnet trainer with extra lace hole allows topline to be adjusted.

Shown in white with black, all black and beige leather with orange.



Airnet colours are available in: orange, blue, grey, black and white



Sprinter

Trainer shoe. Can also be made with antiscuff mudguard & o/s counter for added durability.

Shown in black leather.



Pacer

Trainer boot.

Shown in white/black leather and pink/cream leather.



Abby

Extended opening casual shoe with multi-flower embroidery.

Shown in white leather & soft pink leather.



Sammy

Easy opening casual shoe with single flower embroidery

Shown in white leather and natural leather.



Design your own Sammy!
Choose from three embroidery designs.



Python

Extended opening trainer bootee with contrast animal print flash & collar.

Shown in white leather with natural animal print flash and collar.



Cobra

Extended opening trainer bootee with animal print nubuck and contrast leather flash & collar.

Shown in white leather with animal print facing and soft pink leather flash and collar.



Vyper

Extended opening trainer bootee.

Shown in navy nubuck burgundy flash & sage nubuck with pistachio flash.



Insoles

KIDZ STEP 5 Degree Orthotic

Kidzstep heel cups are designed to support all 3 of the arches. Medial, Lateral and Transverse metatarsal arches. Provides motion control of the foot at each phase of the gait cycle.

Features:

- Medium Density EVA
- Reinforced plastic rear shell
- 5 degree Preposted rearfoot
- Easily ground for practitioner adjustment

Available in sizes

19-20, 21-22, 23-24, 25-26, 27-28, 29-30, 31-32, 33-34



NEUROTEC Orthotic

Neurotec is designed to give support under the medial arch and the peroneal notch, which gives stability to the subtalar joint and destabilises inversion and eversion.

The metatarsal head arch stabilises the forefoot and offloads the metatarsal heads in the forefoot. The heel cup stabilises the heel and the toe crest dorsiflexes the toes to improve tone.

Features:

- Each pair is supplied with 4 degree adhesive backed medial posts
- 1mm adhesive backed EVA build up material

Indications for use:

- Mid neurological deficit affecting the foot
- Habitual reflex deformities
- Hypertonic muscle imbalance



Euro	11	12	13	14	15	16	17	18	19	20
Shell length mm	75	81	87	93	100	103	112	118	126	132
Euro	21	22	23	24	25	26	27	28	29	30
Shell length mm	138	144	151	157	163	169	177	183	189	195



School Collection 21 - 40

Andrea

Normal opening desert boot.
Shown in black patent leather lining.



Chorlton

Extended opening boot with
Anti scuff quarters.
Shown in black leather with
Xymid anti scuff quarters.



Raven

Extended opening boot.
Shown in black leather.



Dodger

Normal opening boot. Shown with
optional commando sole units.
Shown in black leather.



Giselle

Normal opening sandal with tear drop
punch. Double strap velcro fasten.
Shown in black leather and black patent.



Sophie

Normal opening double return
velcro fastening.
Shown in black patent.



Holly

Normal opening sandal with flower
detail. Double strap velcro fasten.
Shown in black leather.





Samson

Normal opening. Plain front shoe.
Shown with optional commando sole units.

Shown in black leather.



Wigwam

Shown with optional
commando sole units.

Shown in black leather.



Order any of
these in Kidztuff
spec!



Storm

Bowling shoe style with anti
scuff toe cap & facings.


Shown in black leather with
anti scuff toe piece



Sapphire

Court shoe style sandal.
Plain front. Velcro fasten.

Shown in grey patent.



Add a black embroidered
rose - just order Opal.



Mohawk

Shown with optional
commando sole units.

Shown in black leather.



Hurricane

Normal Opening. Shown with
optional commando sole units.

Shown in black leather.



Kidztuff Collection

All the shoes on these two pages are available as Kidztuff® specification:

- Extra strong extended heel counters
- Full leather linings
- Rigid hardwearing sole
- Stronger, thicker leathers
- Tough, fashionable soling materials
- Increased backheight for stability

Mercury

Extended opening boot with added counter stability and D-ring fastening

Shown in black leather.



Don't want anti scuff?
Just specify on your order.



Galaxy

Extended opening boot with extra outside counter stability, heavy stitch detail and D-ring fastening.

Shown in black patent with anti scuff toe cap.

Saturn

Kidztuff version of our popular Raven extended opening boot.

Shown in Black patent and black leather.



Pluto

Extended opening boot with Anti scuff toecap and heel

Shown in black leather.



Order any of these
in Standard Reeco
spec. Any colour
combinations.



Jess

Desert boot with caramel rubber outsole.

Shown in brown bruise leather.



Star

Extended opening boot with extra outside quarter stability and D-ring fastening.

Shown in black leather with black daisy and patent flash.



Kidz step 5 degree insole - see page 17

Whatever you want... just ask 23

Moon

2 Strap velcro sandal with padded collar.

Shown in black patent.



Sun

Classic brogue full vamp shoe with leather lining.

Shown in caramel patent and black leather.



Orion

Plain front T-bar shoe

Shown in black leather.



Want a trim?
Add flowers or embroidery.



Earth

Classic brogue with nubuck collar and vamp.

Shown in black leather.

Alex

Fully opening velcro fastening shoe, ideal for AFO or swollen feet.

Shown in navy nubuck.



AFO Collection

Choose any of our Extended opening styles and ask for AFO last.
Here is a small selection...

Raven

Extended opening boot.

Shown in black leather and brown bruiser.



Moon

2 Strap velcro sandal with padded collar.

Shown in brown bruiser leather.

Vyper

Extended opening trainer bootee.

Shown in navy nubuck burgundy flash & sage nubuck with pistachio flash.



Cougar

Extended opening rugged boot with leather contrast collar. Commando soles as standard.

Shown in honey nubuck.

Sammy

Easy opening casual shoe with single flower embroidery

Shown in white leather.



Design your own Sammy!
Choose from three embroidery designs.



Abby

Extended opening casual shoe with multi-flower embroidery.

Shown in white leather & soft pink leather.



Straight & Reverse

All our styles are suitable for straight and reverse last. Below is a small selection...

Ruby

Court shoe style sandal with heart detail. Velcro fasten.

Shown in black leather.



Sapphire

Court shoe style sandal. Plain front. Velcro fasten.

Shown in red patent.



Raven

Extended opening boot.

Shown in black leather.



Vyper

Extended opening trainer bootee.

Shown in sage nubuck with pistachio flash.



Stability & Toe Walkers



Anti-Toe Walker Boots

The ATW Boots are specifically designed to aid those children whose gait develops "toe walking" characteristics.

- Specific Padding in upper, insole and tongue to improve proprioceptive input
- Increased upper height to aid ankle control
- Rigid sole unit with mild rocker
- Commando sole option

Stability Boots

The REED Medical Stability Boot incorporates uniquely designed extended stiffeners and padding which stabilizes the lower talus whilst maintaining a firm yet cushioned fit around the upper talus and malleoli.

- Unique Extended Stiffener design to aid control
- Specific Padding to maximise comfort around the malleoli
- Increased upper height for maximum stability
- Padded tongue
- Commando sole option

The ATW and Stability Boots are both available in the following 4 styles...

Reggae



Disco



Swing



Jazz




Ponseti Boots

REED Medical Ponseti Boots are available for use with adjustable Dennis Browne Bar, The Horton Click Bar and the innovative Dobbs Dynamic Bar.



Topline of the boot is angled as indicated by the Ponseti design



Trimmable Straps



Heel check aperture



Adjustable straps ensure optimum heel retention whilst soft padding to the tongue distributes pressure evenly



Soft insoles provide extra comfort

REED Medical Ponseti Boots are supplied in pairs and are available in the following sizes:

70mm | 80mm | 90mm | 100mm | 110mm | 120mm | 130mm | 140mm | 150mm | 160mm



Dobbs Bars




Allows for easy adjustment of external rotation

The Original Dobbs Bar

The Dobbs bar brace allows the child to move and be active preserves muscle strength in the foot and ankle and is less restrictive than a traditional brace.

Attaches to Markell Shoes / Custom AFOs and Ponseti boots with adaptor (comes with brace when ordered with boots)

- Available in sizes 7" – 9" 9" – 14"
- Adjustable 8"-14" Bar Size.



The Dobbs Bar allows for independent movement and greater mobility

The Dobbs Bar is adjustable from 8 to 14 inches

New Quick Release Spring Assist Dobbs Bar

Quick release mechanism allows Ponseti & AFO's to slide on and off easily & let's children kick and move their legs independently while maintaining abduction and brings the child's feet back to neutral in a resting position.

The spring assist provides a dynamic stretch on the calf allowing patients to build calf muscle strength while moving and playing in the brace which helps mitigate the calf atrophy commonly seen in clubfoot patients.

A spring reset brings the child's feet back to neutral in a resting position. The spring maintains an active dorsiflexion stretch while still allowing the child to play and crawl in the brace.

- Attaches to Ponseti / Mitchell Boots
- Available in sizes 7" – 9" 9" – 14"
- Adjustable 8"-14" Bar Size.



Mitchell Quick Release Spring Assist Dobbs Bar

Quick Release mechanism allows shoes and custom AFO's to snap on and off for changing, carseats, and high chairs.

A spring reset brings the child's feet back to neutral in a resting position. The spring maintains an active dorsiflexion stretch while still allowing the child to play and crawl and build calf muscle strength.

- Attaches to Markell Shoes and Custom AFOs
- Available in sizes 7" – 9" 9" – 14"
- Adjustable 8"-14" Bar Size.




Whatever you want... just ask

Colour chart

Leather



Patent




Nubuck



Suede



Satin



Pearl



Embroidery options

Choose your own embroidery design to be placed on the lateral quarter of your shoe or boot, or on the vamp of the Opal style in place of the single flower.



Ballerina



Heart



Crown



Flamingo



Dragonfly



Multi flower



Single flower



Unicorn



Snake



Dinosaur



Teddy

Accommodation chart

Kids Extra Extra Narrow

Size	Joint Width	Joint Circ.	Instep Circ.
15	43mm	115mm	125mm
16	44mm	120mm	130mm
17	46mm	125mm	135mm
18	47mm	130mm	140mm
19	49mm	135mm	145mm
20	52mm	140mm	150mm
21	53mm	145mm	155mm
22	55mm	150mm	160mm
23	56mm	155mm	165mm
24	58mm	160mm	170mm
25	59mm	165mm	175mm
26	60mm	170mm	180mm
27	62mm	175mm	185mm
28	63mm	180mm	190mm
29	65mm	185mm	195mm
30	67mm	190mm	200mm
31	69mm	195mm	205mm
32	70mm	200mm	210mm
33	71mm	205mm	215mm
34	73mm	210mm	220mm
35	75mm	215mm	225mm
36	76mm	220mm	230mm
37	78mm	225mm	235mm
38	80mm	230mm	240mm
39	82mm	235mm	245mm
40	84mm	240mm	250mm

Kids Standard

Size	Joint Width	Joint Circ.	Instep Circ.
15	54mm	135mm	145mm
16	56mm	140mm	150mm
17	57mm	145mm	155mm
18	59mm	150mm	160mm
19	60mm	155mm	165mm
20	62mm	160mm	170mm
21	63mm	165mm	175mm
22	65mm	170mm	180mm
23	66mm	175mm	185mm
24	68mm	180mm	190mm
25	69mm	185mm	195mm
26	71mm	190mm	200mm
27	72mm	195mm	205mm
28	74mm	200mm	210mm
29	75mm	205mm	215mm
30	77mm	210mm	220mm
31	78mm	215mm	225mm
32	80mm	220mm	230mm
33	81mm	225mm	235mm
34	83mm	230mm	240mm
35	84mm	235mm	245mm
36	87mm	240mm	250mm
37	89mm	245mm	255mm
38	91mm	250mm	260mm
39	92mm	255mm	265mm
40	94mm	260mm	270mm

Kids Extra Narrow

Size	Joint Width	Joint Circ.	Instep Circ.
15	47mm	120mm	130mm
16	48mm	125mm	135mm
17	50mm	130mm	140mm
18	51mm	135mm	145mm
19	53mm	140mm	150mm
20	56mm	145mm	155mm
21	57mm	150mm	160mm
22	59mm	155mm	165mm
23	60mm	160mm	170mm
24	62mm	165mm	175mm
25	63mm	170mm	180mm
26	64mm	175mm	185mm
27	66mm	180mm	190mm
28	67mm	185mm	195mm
29	69mm	190mm	200mm
30	71mm	195mm	205mm
31	73mm	200mm	210mm
32	74mm	205mm	215mm
33	75mm	210mm	220mm
34	77mm	215mm	225mm
35	79mm	220mm	230mm
36	80mm	225mm	235mm
37	82mm	230mm	240mm
38	84mm	235mm	245mm
39	86mm	240mm	250mm
40	88mm	245mm	255mm

Kids Wide

Size	Joint Width	Joint Circ.	Instep Circ.
15	56mm	140mm	150mm
16	58mm	145mm	155mm
17	59mm	150mm	160mm
18	61mm	155mm	165mm
19	63mm	160mm	170mm
20	65mm	165mm	175mm
21	67mm	170mm	180mm
22	69mm	175mm	185mm
23	71mm	180mm	190mm
24	73mm	185mm	195mm
25	74mm	190mm	200mm
26	75mm	195mm	205mm
27	76mm	200mm	210mm
28	78mm	205mm	215mm
29	79mm	210mm	220mm
30	81mm	215mm	225mm
31	83mm	220mm	230mm
32	85mm	225mm	235mm
33	87mm	230mm	240mm
34	89mm	235mm	245mm
35	90mm	240mm	250mm
36	92mm	245mm	255mm
37	93mm	250mm	260mm
38	95mm	255mm	265mm
39	96mm	260mm	270mm
40	98mm	265mm	275mm

Kids Narrow

Size	Joint Width	Joint Circ.	Instep Circ.
15	51mm	130mm	140mm
16	52mm	135mm	145mm
17	54mm	140mm	150mm
18	55mm	145mm	155mm
19	58mm	150mm	160mm
20	60mm	155mm	165mm
21	61mm	160mm	170mm
22	63mm	165mm	175mm
23	64mm	170mm	180mm
24	66mm	175mm	185mm
25	67mm	180mm	190mm
26	68mm	185mm	195mm
27	70mm	190mm	200mm
28	71mm	195mm	205mm
29	73mm	200mm	210mm
30	75mm	205mm	215mm
31	77mm	210mm	220mm
32	78mm	215mm	225mm
33	79mm	220mm	230mm
34	81mm	225mm	235mm
35	83mm	230mm	240mm
36	84mm	235mm	245mm
37	86mm	240mm	250mm
38	88mm	245mm	255mm
39	90mm	250mm	260mm
40	92mm	255mm	265mm

Kids AFO

Size	Joint Width	Joint Circ.	Instep Circ.
15	56mm	145mm	160mm
16	58mm	150mm	165mm
17	59mm	155mm	170mm
18	61mm	160mm	175mm
19	63mm	165mm	180mm
20	65mm	170mm	185mm
21	67mm	175mm	190mm
22	69mm	180mm	195mm
23	71mm	185mm	200mm
24	73mm	190mm	205mm
25	74mm	195mm	210mm
26	75mm	200mm	215mm
27	76mm	205mm	220mm
28	78mm	210mm	225mm
29	79mm	215mm	230mm
30	81mm	220mm	235mm
31	83mm	225mm	240mm
32	85mm	230mm	245mm
33	87mm	235mm	250mm
34	89mm	240mm	255mm
35	90mm	245mm	260mm
36	92mm	250mm	265mm
37	93mm	255mm	270mm
38	95mm	260mm	275mm
39	96mm	265mm	280mm
40	98mm	270mm	285mm

REED

STEP INTO YOUR COMFORT ZONE

FOOTWEAR CATALOGUE

Terms & Conditions

- UK despatch is next day service. Next day timed delivery is available but will incur extra cost
- All adaptations are carried out above the feather edge so as not to alter the bottom pattern so that a standard unit will fit. Should bottom pattern alterations be necessary this will then become bespoke and may incur further costs.
- Children's footwear is supplied with 1 x 3mm footbed.
- Standard sole & heels and through wedge soles are made from EVA.
- Colour reproduction – We make every effort to ensure the colour you see in our catalogue and on our website is accurate, however, colours can be different in photographs and shades can vary from batch to batch and over the life of our catalogue.
- Documentation relating to the confirmation of orders should be clearly marked 'confirmation' to avoid duplication. Duplicated orders through originators error will be chargeable.
- Only Reeco stock footwear is returnable without any alteration / adaption.
- Straight/Reverse & Split sizes & Widths are non-returnable.
- A restocking charge of 30% is chargeable on all returned footwear.
- A return cost of £30.00 is charged for exchange of footwear.
- Returnable footwear must be in a resalable condition and received by us within 1 calendar month following despatch for a refund.
- 1st trial fitting is free of charge and subsequent fittings may be charged for. Accurate measurements are essential when placing your order as re-lasting the upper and new uppers etc following trial fitting will be charged for. An invoice will be raised for all orders not returned for completion with 90 days.
- Split sizes, split fittings, straight/reverse, Afo and alterations to lasts and special orders will be charged in accordance with the current price list and are non-returnable.
- We only use return Velcro fastening straps on all our shoes to enable greater adjustability. Lay on velcro fastening is available as an option but is then non- returnable.
- A restocking charge of 30% is chargeable on all returned footwear. £30.00 is charged for exchange of footwear.
- All shelf stock items; Ponseti boots, bars and insoles are returnable for full refund unused and in a saleable condition.
- The company reserves the right to alter specification or withdraw styles without prior notice. Specification changes are only undertaken when product improvements are identified or a component material is no longer available.
- All goods remain the property of Reed Medical Ltd until paid for in full.
- Website – All our catalogues and brochures are available for download from our website **www.reedmedical.co.uk** our website is a live document and we constantly feed it with information regarding all of our new styles, current news items etc. Please visit it regularly or follow us on facebook and you will receive updates/information regularly and as it happens.
- Should you require further information or wish to speak with our Customer care and Technical team please contact us **on 01254 503333**.

Reed Medical Ltd

Thwaites Close
Shadsworth Business Park
Blackburn Lancashire
United Kingdom BB1 2QQ

Tel 01254 503333
Fax 01254 503360
E-mail admin@reedmedical.co.uk

www.reedmedical.co.uk

