

# reed

bespoke footwear

ISSUE 9

FOOTWEAR  
COLLECTION





# WELCOME

## to the 9th edition of the reed Orthopaedic footwear catalogue

**Reed Medical was established in 1960  
manufacturing bespoke, modular, stock and  
diabetic specification footwear for the NHS.**

We manufacture a range of high-quality footwear to exacting standards. You will find lots of new styles and materials including our Oedema range of stretch shoes, plus our new lightweight and flexible ribbed sole and heel unit.

We have also increased our range of 'off the shelf' footwear for the Diabetic, Rheumatic and Swollen feet, most with stretch or soft material uppers. The Podartis range also includes a Biomechanical outsole with timing rocker.

Our workforce consists of highly skilled and experienced personnel, combining traditional skills with modern production methods and technology to provide functional and stylish footwear.

Our customer service team are friendly and knowledgeable and strive to offer all our customers an unrivalled level of support.



**reed**  
bespoke footwear





# LADIES'

|          |        |
|----------|--------|
| Everyday | 6 - 10 |
|----------|--------|

|         |         |
|---------|---------|
| Leisure | 11 - 13 |
|---------|---------|

|      |         |
|------|---------|
| Boot | 14 - 16 |
|------|---------|

|                |         |
|----------------|---------|
| Stretch/Oedema | 17 - 20 |
|----------------|---------|

|            |    |
|------------|----|
| House Shoe | 21 |
|------------|----|



**KEIRA**

Shown on our classic last in navy leather, velcro option 1 x 38mm.



**RHIANNON**

Shown on our classic last in alabaster leather, velcro option 1 x 38mm.



**THERAPY**

Shown on our classic last in brown nubuck, velcro option 2 x 25mm.



**CARINA**

Shown on our classic last in taupe nubuck, velcro option 1 x 38mm.



**DELICIA**  
Shown on our classic last in burgundy leather/nubuck, velcro option 1 x 38 mm.



**ROSIE**  
Shown on our classic last in taupe leather, velcro option 1 x 38 mm.



**BETH T-BAR**  
Shown on our classic last in black leather, return velcro strap.



**FULLFLEX**  
Shown on our classic last in brown leather/neoprene, velcro option 1 x 38 mm.



**FLEXTOP**

Shown on our classic last in navy leather/neoprene, velcro option 1 x 38mm.



**HAZEL**

Shown on our classic last in beige leather, lace only.



**SUSAN**

Shown on our nova last in burgundy leather, velcro option 2 x 25mm.



**FLEUR**

Shown on our nova last in beige leather, velcro option 2 x 25mm.



**JANE**  
Shown on our nova last in black leather,  
velcro option 1 x 38 mm.



**BUNIONFLEX**  
Shown on our nova last in black  
leather/neoprene, velcro option 2 x 25mm.



**NOVAFLEX**  
Shown on our nova last in brown  
neoprene/leather, velcro option 1 x 38 mm.



**ELLIE**  
Shown on our nova last in burgundy  
leather, velcro option 2 x 25mm.



**ZOE**

Shown on our nova last in black leather, velcro only.



**VIOLET**

Shown on our nova last in white leather, velcro only.



**JUDY**

Shown on our nova last in alabaster leather, velcro only.



**EVE**

Shown on our nova last in white leather/silver metallic, velcro option 2 x 25mm.



### DENISE

Shown on our classic last in white textured softee/silver flash, velcro option 2 x 25mm.



### HEATHER

Shown on our classic last in purple suede/white airnet, velcro option 2 x 25mm.



### LEXI

Shown on our classic last in beige suede/beige airnet, velcro option 2 x 25mm.



### HELENA

Shown on our classic last in white textured softee/silver metallic, velcro option 2 x 25mm.



### VENUS

Shown on our classic last in cream textured softee/gold metallic, velcro option 2 x 25mm.



### SERENA

Shown on our classic last in white textured softee/pink, velcro option 2 x 25mm.



### BIRD

Shown on our classic last in white leather, velcro option 2 x 25mm.



### BILLIE

Shown on our classic last in royal airnet/suede/reflective flash, velcro option 2 x 25mm.



### EMMA

Shown on our nova last in white leather/beige suede, velcro option 2 x 25mm.



### LUNA

Shown on our nova last in white leather/beige suede, velcro option 2 x 25mm.



### PAULA

Shown on our classic last in white textured softee/baby blue, velcro option 3 x 25mm.



### MARY

Shown on our classic last in navy textured softee/white, velcro option 3 x 25mm.

LADIES' BOOT



**ADIE**

Shown on our classic last in brown  
bruiser, velcro option 3 x 25mm.



**CLAIRE**

Shown on our classic last in mid-brown  
leather/tan nubuck, velcro option 3 x  
25mm.



**SALLY**

Shown on our classic last in black  
leather, velcro option 3 x 25mm.



**LOUISE**

Shown on our classic last in black  
leather, velcro option 3 x 25mm.

LADIES' BOOT



**PENNY**  
Shown on our classic last in navy leather, velcro option 3 x 25mm.



**REBECCA**  
Shown on our classic last in burgundy leather, velcro option 3 x 25mm.



**KIRSTEN**  
Shown on our classic last in honey nubuck, velcro option 3 x 25mm.



**FLEXTOP BOOT**  
Shown on our classic last in navy leather/neoprene, velcro option 3 x 25mm.

LADIES' BOOT



**MEGAN**

Shown on our nova last in navy leather, velcro option 3 x 25mm.



**ERIN**

Shown on our nova last in navy nubuck, velcro option 3 x 25mm.



**ARIEL**

Shown on our classic last in black leather/daisy, velcro option 3 x 25mm.



**CASEY**

Shown on our classic last in black leather/rose embroidery, velcro option 3 x 25mm.



### FULLFLEX BOOT

Shown on our classic last in black leather/neoprene, velcro option 3 x 25mm.



### NEOPRENE BOOT

Shown on our classic last in black neoprene, velcro only.



### DAISY

Shown on our classic last in grey stretch/grey suede/silver, velcro option 2 x 25mm.



### DOLLY

Shown on our nova last in black stretch/black suede & deep stretch collar, velcro only.

All shoes on pages 18 & 19 are heat-mouldable



### ASPEN

Shown on our classic last in black PU/  
PU heat-mouldable, velcro only.



### ELM

Shown on our classic last in black  
leather qtrs/heat-mouldable vamp,  
velcro option 1 x 38 mm.



### YEW

Shown on our classic last in black  
PU/Lycra heat-mouldable, velcro only.



### MAPLE

Shown on our classic last in black  
leather, deep neoprene collar/PU  
heat-mouldable vamp, velcro only.



**OAK**

Shown on our classic last in taupe leather, deep neoprene collar/PU heat-mouldable vamp, velcro only.





**BEECH**

Shown on our classic last in black leather, deep neoprene collar/PU heat-mouldable vamp, velcro only.




**LARCH**

Shown on our classic last in black leather, deep neoprene collar/PU heat-mouldable vamp, velcro only.



Unique heat-mouldable vamp for localised customisation e.g. bunion width or clawed toes, etc. An insert is used to bellow out the material in specific areas while heat is applied with a low temperature heat gun or hair dryer. Once cooled the material will hold its shape until further heat is applied.

**PU HEAT-MOULDABLE**



Black

Navy

Brown

Burgundy

Taupe



### VANGUARD SHOE

Shown on our nova last in black velcro receptive. Black only, velcro only.



### VANGUARD BOOT

Shown on our nova last in black velcro receptive. Black only, velcro only.



### DONCASTER BOOT

Shown on our classic last in black velcro receptive/Lycra heat-mouldable vamp, velcro only.

LADIES' HOUSE SHOE



**PISCES**  
Shown on our classic last in pink felt/fur lining, velcro only.



**GEMINI**  
Shown on our classic last in beige felt/fur lining, velcro only.



**LIBRA**  
Shown on our classic last in black felt/fur lining, velcro only.



**VIRGO**  
Shown on our classic last in navy felt/deep neoprene collar, velcro only.



# MEN'S

Everyday 24 - 27

Leisure 28 - 31

Boot 32 - 34

Stretch/Oedema 35 - 38

House Shoe 39



**CALLUM**

Shown on our classic last in black leather, velcro option 1 x 38mm.



**THOMAS**

Shown on our classic last in brown leather, velcro option 1 x 38mm.



**AUSTEN**

Shown on our classic last in brown nubuck, velcro option 2 x 25mm.



**THERAPY**

Shown on our classic last in brown nubuck, velcro option 2 x 25mm.



**CONNOR**  
Shown on our classic last in mid-brown leather, velcro option 1 x 38mm.



**FULLFLEX**  
Shown on our classic last in brown leather/neoprene, velcro option 1 x 38 mm.



**THOR**  
Shown on our classic last in brown bruiser, velcro option 1 x 38mm.



**LOGAN**  
Shown on our classic last in brown bruiser/red, velcro option 2 x 25mm.



**FLEXTOP**

Shown on our classic last in navy leather/neoprene, velcro option 1 x 38mm.



**TUCSON**

Shown on our nova last in mid-brown leather, velcro option 1 x 38mm.



**CARSON**

Shown on our nova last in brown nubuck, velcro option 1 x 38mm.



**NOVAFLEX**

Shown on our nova last in brown leather/neoprene, velcro option 1 x 38mm.



**BUNIONFLEX**  
Shown on our nova last in black leather/neoprene, velcro option 2 x 25mm.



**MONTANA**  
Shown on our nova last in brown bruiser, velcro option 2 x 25mm.



**RENO**  
Shown on our nova last in brown bruiser contrast stitching, velcro option 1 x 38mm.



### **DALY**

Shown on our classic last in black textured softee/white, velcro option 2 x 25mm.



### **BIRD**

Shown on our classic last in white leather, velcro option 2 x 25mm.



### **BOLT**

Shown on our classic last in white textured softee/red, velcro option 2 x 25mm.



### **COE**

Shown on our classic last in navy textured softee/white, velcro option 2 x 25mm.



### **MURRAY**

Shown on our classic last in navy suede/navy airnet, velcro option 2 x 25mm.



### **CHRISTIE**

Shown on our classic last in navy leather/suede, velcro option 2 x 25mm.



### **BILLIE**

Shown on our classic last in royal airnet/suede/reflective flash, velcro option 2 x 25mm.



### **DASH**

Shown on our classic last in black suede/white, velcro option 2 x 25mm.



### BOTHAM

Shown on our classic last in tan leather, velcro option 2 x 25mm.



### TANGO

Shown on our classic last in red leather/white, velcro option 2x25mm



### LENNOX

Shown on our classic last in navy textured softee/white, velcro option 3 x 25mm.



### BAYLOR

Shown on our classic last in black leather, velcro option 3 x 25mm.



**LUKE**  
Shown on our classic last in aqua green airnet/aqua green suede, velcro option 3 x 25mm.



**EMBREY**  
Shown on our classic last in black leather/grey suede, velcro option 3 x 25mm.



**BROADY**  
Shown on our classic last in navy airnet/suede, velcro option 3 x 25mm.



**OVETT**  
Shown on our classic last in white textured softee/navy, velcro option 3 x 25mm.

## MEN'S BOOT



### BEVAN

Shown on our classic last in brown leather, velcro option 3 x 25mm.



### FLEXTOP BOOT

Shown on our classic last in navy leather/neoprene, velcro option 3 x 25mm.



### TYSON

Shown on our classic last in brown bruiser, velcro option 3 x 25mm.



### KYLE

Shown on our classic last in honey nubuck/brown leather, velcro option 3 x 25mm.

## MEN'S BOOT



### LORCAN

Shown on our classic last in black leather, velcro option 3 x 25mm.



### SAXON

Shown on our classic last in black leather, velcro option 3 x 25mm.



### WAVE

Shown on our classic last in brown bruiser/brown leather, velcro option 3 x 25mm.



### WADE

Shown on our classic last in brown bruiser/brown leather, velcro option 3 x 25mm.

## MEN'S BOOT



### GABLE

Shown on our classic last in black leather/anti-stuff toe/neoprene collar, velcro option 3 x 25mm.



### JOSEPH

Shown on our nova last in navy suede, velcro option 2 x 25mm.



### MARK

Shown on our nova last in brown bruiser, velcro option 2 x 25mm.



### NEVADA

Shown on our nova last in mid-brown leather, velcro option 3 x 25mm.



### FULLFLEX BOOT

Shown on our classic last in black leather/neoprene, velcro option 3 x 25mm.



### NEOPRENE BOOT

Shown on our classic last in black neoprene, velcro only.



### FULLFLEX

Shown on our classic last in brown leather/neoprene, velcro option 1 x 38mm.



### BLAKE

Shown on our classic last in black leather/neoprene, lace option available.

All shoes on pages 36 & 37 are heat-mouldable



### ELM

Shown on our classic last in brown leather/PU heat-mouldable, velcro option 1 x 38mm.



### YEW

Shown on our classic last in black leather/Lycra heat-mouldable, velcro only.



### MAPLE

Shown on our classic last in black leather/PU heat-mouldable vamp, velcro only.



### OAK

Shown on our classic last in black leather, neoprene/PU heat-mouldable, velcro only.





**ASPEN**  
Shown on our classic last in black  
PU/PU heat-mouldable, velcro only.



**LARCH**  
Shown on our classic last in black  
leather/PU heat-mouldable, velcro only.




**BEECH**  
Shown on our classic last in black  
leather, deep neoprene collar/PU heat-  
mouldable, velcro only.



Unique heat-mouldable vamp for  
localised customisation e.g. bunion  
width or clawed toes, etc. An insert  
is used to bellow out the material in  
specific areas while heat is applied with  
a low temperature heat gun or hair dryer.  
Once cooled the material will hold its  
shape until further heat is applied.

**PU HEAT MOULDBALE**





### VANGUARD SHOE

Shown on our nova last in black velcro receptive. Black only, velcro only.



### VANGUARD BOOT

Shown on our nova last in black velcro receptive. Black only, velcro only.



### DONCASTER BOOT

Shown on our classic last in black velcro receptive/Lycra heat-mouldable vamp, velcro only.



**TAURUS**  
Shown on our classic last in black felt/fur lining, velcro only.



**LEO**  
Shown on our classic last in black felt/fur lining, velcro only.



**VIRGO**  
Shown on our classic last in navy felt & neoprene, velcro only.



# OFF THE SHELF

## PODARTIS

|               |    |
|---------------|----|
| Deambulo Line | 42 |
|---------------|----|

|               |         |
|---------------|---------|
| Activity Line | 43 - 45 |
|---------------|---------|

|                          |         |
|--------------------------|---------|
| Diabetic & Classic shoes | 46 - 47 |
|--------------------------|---------|

|                         |         |
|-------------------------|---------|
| Acute Phase Off-Loading | 48 - 49 |
|-------------------------|---------|

|          |         |
|----------|---------|
| Podowell | 50 - 51 |
|----------|---------|

|                   |    |
|-------------------|----|
| Lymphoedema Shoes | 52 |
|-------------------|----|



## DEAMBULO L

SA200 Fitting 12 Variable Size 39-47  
Stretchable upper in Setaflex. Lightweight anti-stumble flexible outsole.



## DEAMBULO X

SA400 Fitting Variable Size 36-47  
Fully opening with variable volume. Biomechanical anti-stumbling outsole.



## DEAMBULO SPORT

SA300 Fitting 12 Variable Size 35-41  
Stretchable upper over forefoot and dorsum in Setaflex. Ultra lightweight anti-stumble flexible outsole.



## FLY

SA571 Fitting 12 Variable Size 35-41  
Flexpell® stretchable, breathable upper and patent leather with elastic strap. Flexible anti-stumble biomechanical outsole 20mm heel height.



**FANCY SILVER**

PA20178 Fitting 12 Size 35-41  
Soft Flexpell® upper over toes.  
Biomechanical semi-rigid EVA Vibram®  
outsole. Timing Rocker® sole.



**FANCY BLACK**

PA20171 Fitting 12 Size 35-41  
Soft Flexpell® upper over toes.  
Biomechanical semi-rigid EVA Vibram®  
outsole. Timing Rocker® sole.



**EVO MIDNIGHT VELCRO**

PA20418 Fitting 12.5 Size 35-41  
Super soft Flexpell® upper over toes.  
Biomechanical semi-rigid Piuma® high  
abrasion resistance. Timing Rocker® sole.



**ADELLE SILVER VELCRO**

SR713118 Fitting 12.5 Size 35-41  
Super soft Flexpell® upper over toes.  
Biomechanical anti-stumbling outsole.  
Timing Rocker® sole.



## FREEDOM EVO

PA34508 Unisex Fitting 12 Size 35-47  
Water resistant nappa stretch.  
Semi-rigid biomechanical Timing Rocker® sole Piuma®.



## ACTIVITY ONE EVO VELCRO

PA30501 Fitting 12.5 Size 39-47  
Wide forepart for deformities,  
swelling or diabetes. Wellwalk® energy  
saver technology outsole.



## HERO EVO

PA30030 Fitting 12 Size 39-46  
Stretchable upper to accommodate  
deformities. Biomechanical anti-stumbling  
Timing Rocker® Piuma® outsole.

## DCS TECHNOLOGY

The DCS (Dynamic Cross System®) is a composite fibre insert located inside the Vibram® tread that guarantees:

- Controlling excess of pronation and supination (tested up to 140kg)
- Increased lift to avoid sole collapse in overweight and obese people
- Elastic smooth transition: very high elastic response to facilitate the unfolding of the podalic helix and reduce energy demand.

**Clinical Study.** Evaluation of the effectiveness of Activity DCS-AFO in patients with stroke outcomes.

Rehabilitation Centre 'villa beretta', Costa Masnaga (Dr. F. Molteni)



### WILD

PA365521 Fitting 14.5 Size 39-47  
Rigid biomechanical Timing Rocker® sole & extended stabilising heel counter.



### SPORT BLACK

PA37781 Fitting 13 Size 39-47  
Resistance against pronation and supination uses DCS (Dynamic Cross System®) technology.



### ITALIA BLACK

PA4131 Fitting 12 Size 35-41  
Semi-rigid Timing Rocker® biomechanical outsole uses DCS (Dynamic Cross System®) technology. Carbon fibre insert to curb excessive pronation/supination.



### GIOTTO

SD3001 Fitting 11.5 Size 39-47  
Semi-rigid biomechanical EVA sole with Vibram® outsole. Heat-mouldable upper.



### BOTERO 11 VELCRO

SD370291 Fitting 11.5 Size 39-46  
Flexpell® stretchable upper to accommodate forefoot. Semi-rigid biomechanical EVA sole with Vibram® outsole.



### LADY ALFA

SD690191 Fitting 14 Size 35-41  
Leather upper with Flexpell® stretchable vamp. Biomechanical rigid EVA outsole.



### GIOTTO GV

SD432101 Fitting 14 Size 39-46  
Heat-mouldable upper. Semi-rigid biomechanical outsole powered by energy saver technology. Heel strike 20%, Midstance 55%, Toe off 25%.



**X-DIAB 14**

SD30091 Fitting 14 Size 39-46  
Leather upper with Flexpell® stretchable vamp. Wellwalk® biomechanical outsole and Vibram® energy saver technology.



**X-DIAB 16**

SD3471 Fitting 16 Size 36-46  
Leather upper with Flexpell® stretchable vamp. Wellwalk® biomechanical outsole and Vibram® energy saver technology.



**Since 2017, Podiapp has been the reference application for the care and treatment of the sensitive foot.**

**Now completely renewed, with many new functions and tips for the well-being of all feet!**

Scan me!



[www.reedmedical.co.uk](http://www.reedmedical.co.uk)



Stabil D + Modus insole  
Lab Novel, 2007

## STABIL D


PO700 Size 35-48  
Removable insole to insert MODUS or custom-made insole. Rigid carbon fibre insole, rigid outsole and lateral stabilising removable whale bones.



## TD-FOREFOOT & MIDFOOT


PO500 Size 35-48  
No crutch: the biomechanical outsole enables the patient to walk safely and get back to activity sooner. Removable insole to enable use with bespoke or our MODUS insole for modular off-loading.

# ACUTE PHASE OFF-LOADING



## WPS-FOREFOOT OFF-LOADER

PO100 Size 35-48  
Stability: ideal angle of ankle dorsiflexion to stabilise the walk and reduce “step effect”. Rigid counters and rear stabilisers to ensure maximum stability. Removable EVA insole allows for inclusion of MODUS or custom-made insole.



## TERAHEEL® HEEL PRESSURE OFF-LOADING

PO400 Size 35-46  
Full opening with velcro fastening. Side counters. Off-loading insole in soft PU.



## FOR USE WITH ALL OFF-LOADING PRODUCTS

### MODUS INSOLE

PL345 Size 35-48  
Covered in dermocompatible DIAPOD, antibacterial and antifungal properties. Ideal for plantar surface ulcer pressure relief.

**Podowell® offers a range of roomy, flexible and washable footwear.**



SANITIZED® Treatment that limits the development of bacteria and odours even after washing



Hand-wash at 30°C maximum



Stretchable material on sensitive zones to protect sore feet



**Extra Extra Large**  
Maximum of volume to fit oedema or lymphoedema.



**ALTITUDE**

Size 36-48  
Ideal for use with oedematous or neuropathic feet, ulceration, incontinent users or swollen feet, and will accommodate bandaged feet.



**XAVIER**

Size 36-48  
Ideal for use with oedematous or neuropathic feet, ulceration, incontinent users or swollen feet, and will accommodate bandaged feet.



**Extra Large**

Substantial fitting room to fit post-operative, diabetic or seriously swollen feet.



**ALADIN WIDE FIT**

Size 36-48

A washable closed toe soft stretch vamp shoe, extra roomy with removable inlay.

**ALIX WIDE FIT**

Size 35-46

Microflex vamp, hand-washable, wide fitting shoe with removable insole.



**ANDALOU WIDE FIT**

Size 35-46

Open toe, soft textile upper with antibacterial finish and removable inlay for use with a captive Total Contact insole. 6mm removable insole to replace with own.

**ATHOS WIDE FIT**

Size 36-48

A closed toe soft shoe, extra roomy with removable inlay.

LYMPHOEDEMA SHOES



**LYMPH VELCRO 60330**  
Designed to accommodate bandaged legs and feet with a wide tongue, a lightweight sole and cushioned insole. The vamp/tongue can be removed to produce a closed sandal.



**LYMPH LACE 60310**  
Designed to accommodate bandaged legs and feet with a wide tongue, a lightweight sole and cushioned insole.



**LYMPH SANDAL 60320**  
The Lymph sandal is made with a stretchable material which fully opens to allow for ease of donning. Includes removable toe/tongue. Designed to accommodate bandaged legs and feet with a wide tongue, a lightweight sole and cushioned insole.

| SIZE CHART FOR LYPH SHOE/SANDAL ONLY |      |     |      |     |      |     |      |     |      |     |      |     |      |      |
|--------------------------------------|------|-----|------|-----|------|-----|------|-----|------|-----|------|-----|------|------|
| SIZE                                 | 36   | 37  | 38   | 39  | 40   | 41  | 42   | 43  | 44   | 45  | 46   | 47  | 48   | 50   |
| FOOT LENGTH                          | 239  | 246 | 253  | 259 | 266  | 273 | 279  | 286 | 293  | 299 | 306  | 313 | 319  | 333  |
| JOINT WIDTH                          | 91   | 93  | 95   | 97  | 99   | 101 | 103  | 105 | 107  | 109 | 111  | 113 | 115  | 119  |
| HEEL WIDTH                           | 60.5 | 62  | 63.5 | 65  | 66.5 | 68  | 69.5 | 71  | 72.5 | 74  | 75.5 | 77  | 78.5 | 81.5 |
| TOE DEPTH                            | 28.5 | 29  | 29.5 | 30  | 30.5 | 31  | 31.5 | 32  | 32.5 | 33  | 33.5 | 34  | 34.5 | 35.5 |
| JOINT CIRC.                          | 270  | 275 | 279  | 284 | 288  | 293 | 297  | 302 | 307  | 311 | 315  | 319 | 323  | 331  |
| INSTEP CIRC.                         | 280  | 285 | 289  | 294 | 298  | 303 | 307  | 312 | 317  | 321 | 325  | 329 | 333  | 341  |



# CUSTOM PRODUCTS

Custom Orthotics 54

Orthotic Prescription Form 55

Our Bespoke Range 56

Casting for Custom Footwear 57


Other Products 58

## REED CUSTOM INSOLES

Technology has always been an integral part of our business at Reed Medical, whether it is last carving or upper design we have always strived to create the perfect blend of both traditional and modern methods to create high-quality products.

It is with great excitement that we continue this tradition by announcing the CAD CAM insole range '**REED CUSTOM**'. Available in Multi Density EVA, Polyurethane, Polypropylene and Carbon Fibre with multi choice coverings.

Once your orthotics are milled, they are expertly hand-finished by a skilled technician to your exact requirements, ensuring quality and perfect fit every time.



# ORTHOTIC PRESCRIPTION FORM

Name:  Date:   
Hospital/Clinic:  Require Date:   
NHS Order Ref:  Clinician:

Make to REEDS Footwear: ☐ Other Footwear: ☐ Size:   
Size:  Width:  Template Supplied: ☐  
Please Draw On Reverse  
Footwear Order Number:   
Full Contact: ☐ Standard: ☐ Pair: ☐ Left Only: ☐ Right Only: ☐  
Unrectified IMP Box/Cast 10% Arch Reduction  
Sulcus Flattening Only  $\frac{3}{4}$  Flattening

## Base Materials:

**EVA:** LD A30 ☐ MD A40 ☐ MD A50 ☐ HD A60 ☐  
PU ☐ Carbon ☐ Polypropylene ☐ Dual ☐

## Base Thickness:

3mm ☐ 6mm ☐ Other:  mm

## Length:

$\frac{3}{4}$  ☐ Sulcus ☐ Full ☐

## Flexibility/Arch Fill:

Flexible ☐ Semi Flex ☐ Rigid ☐

**Posting:** Intrinsic ☐ Extrinsic ☐

**Rearfoot:** Degrees:

**Left -** Medial ☐ Lateral ☐

**Right -** Medial ☐ Lateral ☐

**Forefoot:** Degrees:

**Left -** Medial ☐ Lateral ☐

**Right -** Medial ☐ Lateral ☐

## Off-Loading:

Cut Out ☐ Sink ☐  
Through Base Top

Fill:   
Specify Material

Unfilled: ☐

Please draw on images or  
mark on impression box

## Additional Information:



**Heel Raise:** Short ☐ Long ☐

Left:  mm

Right:  mm

## Kirby Skive:

Degrees:

**Left -** Medial ☐ Lateral ☐

**Right -** Medial ☐ Lateral ☐

## Coverings:

Middle:

Specify Material:

Thickness: 1mm ☐ 3mm ☐ 6mm ☐

Other:  N/A ☐

Top:

Specify Material:

Thickness: 1mm ☐ 3mm ☐ 6mm ☐

Other:  N/A ☐

# OUR BESPOKE RANGE

We provide fully handmade shoes to your cast and draft: see casting for custom footwear opposite then choose a suitable style from this catalogue.

Should you require further information at this point please contact us **01254 503333**



The cast and draft is an important aspect of the custom footwear process. A quality cast and draft is essential for ensuring the fit of the footwear.

An accurate cast will capture the position of the foot/leg and any oedema, and any other orthotic devices that will be worn inside the shoe. The position of the foot/leg is important especially when accommodating a fixed deformity. If the foot is in equinus or plantar flexion, the cast should be taken in the position you require the finished footwear.

### CASTING POSITION

- 1 Whenever possible, cast the patient sitting with their feet flat on the ground and the lower leg perpendicular (at 90 degrees) to the floor.
- 2 For patients with severe rigid deformities, cast them in the position of function. Give the shoemaker as much information as possible on the prescription form, and mark alignment lines on the cast. It is important to ensure the cast reflects the shape you require for your patient. **Correcting the leg alignment once the cast is filled is not the responsibility of the shoemaker, having to rectify the cast is extremely time-consuming.**
- 3 The goal is to create a semi-weight-bearing cast that captures the full width of the metatarsals; the full length of the extended foot and toes; and the contours of the foot at the heel and arch. If the patient is going to wear an AFO, then cast them with it on or send with the cast for the shoemaker to make allowance inside the footwear. If you are unable to cast with AFO on then take your measure over the AFO whilst being worn by the patient and draw around the foot of the AFO. **If an AFO is not available at the time of casting, it is best to send it on to us so that we can incorporate it into the last/cast.**  
  
If a patient is to be cast in a wheelchair or lying down or casts that are taken non-weight-bearing then this should be clearly noted and the angle of leg clearly marked on the cast.
- 4 When pulling on or using the casting stockings, make sure that the toes remain fully extended; otherwise, the footwear may come back too short. It is important to ensure clear toe definition.
- 5 Maintaining the required angle, mark the alignment lines, make pencil marks across the cut-line of the cast to indicate alignment before removing or opening with a cast saw. This is very helpful for the shoemaker when sealing to make a positive mould. Peel the cast away gently so as not to distort the shape.
- 6 Record purchase order number or patient ID and any areas of concern on the cast. After removing, check the final cast and ensure that any areas of concern inside the cast are marked.



After casting is complete, and it is possible, take a draft and measures and/or foot impression. If the patient has an existing comfortable pair of custom footwear, trace the insoles.

A good cast/draft will ensure a pair of well-fitting and comfortable custom footwear.

Examples of sole adaptations available



Through Raise



Socket



Buttress



Custom-made BK Irons



Bespoke Chair boots

- Leather upper
- Sheepskin lining
- Sewn in sock
- Velcro fasten



3/4 Gaitor Arch Support

- Cushioned heel
- Plantar fasciitis
- Heel spurs
- Fallen arches
- Shin splints



Pedag Classic Insoles

- Leather covered soft foam insole
- Arch support
- Heel cushion
- Met button



# TECHNICAL INFORMATION

|                |    |
|----------------|----|
| Last Adaptions | 60 |
|----------------|----|

|                             |    |
|-----------------------------|----|
| Measuring for Reed Footwear | 61 |
|-----------------------------|----|

|                    |    |
|--------------------|----|
| Ladie's Size Chart | 62 |
|--------------------|----|

|                  |    |
|------------------|----|
| Men's Size Chart | 63 |
|------------------|----|

|           |    |
|-----------|----|
| Materials | 64 |
|-----------|----|


|                |    |
|----------------|----|
| Soling Options | 65 |
|----------------|----|

|       |    |
|-------|----|
| Ts&Cs | 66 |
|-------|----|

|                |    |
|----------------|----|
| Kids' Footwear | 67 |
|----------------|----|

## MODULAR LAST ADAPPTIONS

- 1 Toe depth
- 2 Joint width lateral/medial
- 3 Joint depth
- 4 Instep depth
- 5 Heel width (around the heel ball)
- 6 Through depth



Additions are available in multiples of 3mm. Any adaption over 6mm will require a bespoke pattern.

## LAST SHAPES

All of our lasts are pitched for a 25mm heel height.



Our original **Classic** last is an almond toe shape with a rounded toe profile and through depth for 3x3mm inlay. This last is available with **EVA Pyramid** or **PU Ribbed** sole unit.

Our **Nova** last has a modern shaped toe with a square walled toe profile and through depth for 3x3mm inlay. This last is available with our **Hidden depth** PU sole unit.

Our **Nature** last is ideal for the 'at risk' foot. This last has a wider forepart and complements the natural contours of the foot and through depth for 3x3mm inlay. This last is available with **EVA Pyramid** or **PU Ribbed** sole unit.

**Lowline** is available across all above last shapes which has a through depth of 1x3mm inlay.

**Inlays are removable to allow for fitting adjustment.**

**EVA is used for Standard sole and heels, through wedge soles and soles requiring raises, sockets, etc.**

## SHOE SIZE

Measure the length of the foot using your Reed measuring stick. Place the fixed jaw of the stick at the heel and measure to the longest toe. Hold the sliding jaw in place and read the corresponding size on the stick. Add 1½ sizes if patient standing or 2 if sitting. This is the shoe size.

## JOINT WIDTH

Open the jaw on your stick place either side of the joint, hold the jaw in place lift away from the foot and read the measure.

## TOE DEPTH

Place fixed jaw under toes and measure depth to highest toe, hold the jaw in place lift away from the foot and read the measure.

## JOINT CIRCUMFERENCE

Place your tape measure under the met heads and measure around the widest part of the foot.

## INSTEP CIRCUMFERENCE

Place your tape measure under the instep at the waist, meeting at the navicular.

## LEG CIRCUMFERENCE - A MEASURE (BOOTS ONLY)

Place your tape measure at 130mm height and measure around the leg or at the height you wish the boot to be.

## ANKLE CIRCUMFERENCE - B MEASURE (BOOTS ONLY)

Place your tape measure at 90mm height and measure around the ankle.

## INSOLE ALLOWANCE

Consider the depth of insole to be used. All footwear has 9mm allowance excluding Lowline last which has an insole allowance of 3mm.

## FITTING CHARTS

Refer to fitting chart to find the closest match for your measures. **We can't reduce a last.** If you require increased measures, specify your modular adaptations.

All of our footwear is supplied with 3x3mm inlays excluding our Lowline last which has 1x3mm inlay. You will find 2 measures on our accommodation charts e.g. 230 – 243; please note that the higher figure is **WITHOUT** insole allowance.

The footwear is supplied with sole as catalogue; if you require an alternative, please specify when you order.

Standard leg measures on boots will be supplied unless custom measurements are specified. When supplying these please specify the height the measure was taken.

**A:** Circumference of leg at 130mm  
**B:** Circumference at ankle bone

### A&B Measurement example


| Left    | Right |
|---------|-------|
| 241mm A | 228mm |
| 266mm B | 254mm |

If wearing with an AFO please measure with patient wearing and send draft of bottom plate for length.

We offer a free trial fitting on all of our footwear, allowing you to check your fit. Alterations following fitting will incur added costs.

All our footwear is made to order, therefore only our 'off the shelf' products are returnable unworn and box intact: please do not deface the packaging with patient name or ID.

## LAST MEASUREMENT POSITIONS



# LADIES' ACCOMODATION CHART

## LADIES' NARROW LOWLINE

3mm Insole allowance 1x3mm footbed

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 1    | 77mm        | 200mm - 214mm       | 190mm - 200mm        | 28mm           | 54mm       |
| 2    | 79mm        | 206mm - 220mm       | 196mm - 206mm        | 29mm           | 55mm       |
| 3    | 81mm        | 211mm - 226mm       | 200mm - 212mm        | 30mm           | 56mm       |
| 4    | 83mm        | 217mm - 232mm       | 2406mm - 220mm       | 31mm           | 57mm       |
| 5    | 85mm        | 224mm - 238mm       | 212mm - 227mm        | 32mm           | 58mm       |
| 6    | 87mm        | 230mm - 244mm       | 220mm - 232mm        | 33mm           | 59mm       |
| 7    | 89mm        | 236mm - 250mm       | 227mm - 239mm        | 34mm           | 60mm       |
| 8    | 91mm        | 242mm - 256mm       | 232mm - 246mm        | 35mm           | 61mm       |
| 9    | 93mm        | 250mm - 262mm       | 241mm - 253mm        | 36mm           | 62mm       |

## LADIES' NARROW

9mm Insole allowance 3x3mm footbed

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 1    | 77mm        | 209mm - 220mm       | 204mm - 214mm        | 31mm           | 54mm       |
| 2    | 79mm        | 215mm - 226mm       | 210mm - 220mm        | 32mm           | 55mm       |
| 3    | 81mm        | 221mm - 232mm       | 214mm - 226mm        | 33mm           | 56mm       |
| 4    | 83mm        | 227mm - 238mm       | 218mm - 232mm        | 34mm           | 57mm       |
| 5    | 85mm        | 233mm - 244mm       | 224mm - 238mm        | 35mm           | 58mm       |
| 6    | 87mm        | 239mm - 252mm       | 232mm - 244mm        | 36mm           | 59mm       |
| 7    | 89mm        | 245mm - 258mm       | 238mm - 250mm        | 37mm           | 60mm       |
| 8    | 91mm        | 251mm - 264mm       | 242mm - 256mm        | 38mm           | 61mm       |
| 9    | 93mm        | 257mm - 270mm       | 248mm - 262mm        | 39mm           | 62mm       |

## LADIES' STANDARD LOWLINE

3mm Insole allowance 1x3mm footbed

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 1    | 87mm        | 214mm - 224mm       | 202mm - 212mm        | 28mm           | 55mm       |
| 2    | 89mm        | 220mm - 230mm       | 208mm - 218mm        | 29mm           | 56mm       |
| 3    | 91mm        | 226mm - 236mm       | 214mm - 224mm        | 30mm           | 57mm       |
| 4    | 93mm        | 232mm - 242mm       | 220mm - 230mm        | 31mm           | 58mm       |
| 5    | 95mm        | 238mm - 248mm       | 226mm - 236mm        | 32mm           | 59mm       |
| 6    | 97mm        | 244mm - 254mm       | 232mm - 242mm        | 33mm           | 60mm       |
| 7    | 99mm        | 250mm - 260mm       | 238mm - 248mm        | 34mm           | 61mm       |
| 8    | 101mm       | 256mm - 266mm       | 246mm - 256mm        | 35mm           | 62mm       |
| 9    | 103mm       | 262mm - 272mm       | 252mm - 262mm        | 36mm           | 63mm       |

## LADIES' STANDARD

9mm Insole allowance 3x3mm footbed

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 1    | 87mm        | 222mm - 231mm       | 209mm - 222mm        | 31mm           | 55mm       |
| 2    | 89mm        | 228mm - 235mm       | 216mm - 229mm        | 32mm           | 56mm       |
| 3    | 91mm        | 235mm - 248mm       | 222mm - 235mm        | 33mm           | 57mm       |
| 4    | 93mm        | 241mm - 254mm       | 228mm - 241mm        | 34mm           | 58mm       |
| 5    | 95mm        | 248mm - 260mm       | 235mm - 248mm        | 35mm           | 59mm       |
| 6    | 97mm        | 254mm - 266mm       | 241mm - 254mm        | 36mm           | 60mm       |
| 7    | 99mm        | 260mm - 273mm       | 248mm - 260mm        | 37mm           | 61mm       |
| 8    | 101mm       | 266mm - 279mm       | 254mm - 266mm        | 38mm           | 62mm       |
| 9    | 103mm       | 272mm - 286mm       | 260mm - 272mm        | 39mm           | 63mm       |

## LADIES' WIDE LOWLINE

3mm Insole allowance 1x3mm footbed

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 1    | 93mm        | 225mm - 235mm       | 214mm - 224mm        | 28mm           | 56mm       |
| 2    | 95mm        | 231mm - 241mm       | 220mm - 230mm        | 29mm           | 57mm       |
| 3    | 97mm        | 237mm - 247mm       | 226mm - 236mm        | 30mm           | 58mm       |
| 4    | 99mm        | 243mm - 253mm       | 232mm - 242mm        | 31mm           | 59mm       |
| 5    | 100mm       | 249mm - 259mm       | 238mm - 248mm        | 32mm           | 60mm       |
| 6    | 103mm       | 255mm - 265mm       | 244mm - 254mm        | 33mm           | 61mm       |
| 7    | 105mm       | 261mm - 271mm       | 250mm - 260mm        | 34mm           | 62mm       |
| 8    | 107mm       | 267mm - 277mm       | 256mm - 266mm        | 35mm           | 63mm       |
| 9    | 109mm       | 273mm - 283mm       | 262mm - 272mm        | 36mm           | 64mm       |

## LADIES' WIDE

9mm Insole allowance 3x3mm footbed

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 1    | 93mm        | 235mm - 248mm       | 216mm - 228mm        | 31mm           | 56mm       |
| 2    | 95mm        | 241mm - 254mm       | 222mm - 235mm        | 32mm           | 57mm       |
| 3    | 97mm        | 248mm - 260mm       | 228mm - 241mm        | 33mm           | 58mm       |
| 4    | 99mm        | 254mm - 266mm       | 235mm - 248mm        | 34mm           | 59mm       |
| 5    | 100mm       | 260mm - 273mm       | 241mm - 254mm        | 35mm           | 60mm       |
| 6    | 103mm       | 266mm - 279mm       | 248mm - 260mm        | 36mm           | 61mm       |
| 7    | 105mm       | 273mm - 287mm       | 254mm - 266mm        | 37mm           | 62mm       |
| 8    | 107mm       | 279mm - 293mm       | 260mm - 273mm        | 38mm           | 63mm       |
| 9    | 109mm       | 285mm - 300mm       | 267mm - 279mm        | 39mm           | 64mm       |

## LADIES' DEEP

9mm Insole allowance 3x3mm footbed

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 1    | 87mm        | 235mm - 248mm       | 216mm - 228mm        | 40mm           | 55mm       |
| 2    | 89mm        | 241mm - 254mm       | 222mm - 235mm        | 41mm           | 56mm       |
| 3    | 91mm        | 248mm - 260mm       | 228mm - 241mm        | 42mm           | 57mm       |
| 4    | 93mm        | 254mm - 266mm       | 235mm - 248mm        | 43mm           | 58mm       |
| 5    | 95mm        | 260mm - 273mm       | 241mm - 254mm        | 44mm           | 59mm       |
| 6    | 97mm        | 266mm - 279mm       | 248mm - 260mm        | 45mm           | 60mm       |
| 7    | 99mm        | 273mm - 286mm       | 254mm - 266mm        | 46mm           | 61mm       |
| 8    | 101mm       | 279mm - 292mm       | 260mm - 273mm        | 47mm           | 62mm       |
| 9    | 103mm       | 285mm - 298mm       | 266mm - 279mm        | 48mm           | 63mm       |

## LADIES' WIDE & DEEP

9mm Insole allowance 3x3mm footbed

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 1    | 93mm        | 246mm - 260mm       | 234mm - 248mm        | 40mm           | 56mm       |
| 2    | 95mm        | 252mm - 266mm       | 240mm - 254mm        | 41mm           | 57mm       |
| 3    | 97mm        | 258mm - 272mm       | 246mm - 260mm        | 42mm           | 58mm       |
| 4    | 99mm        | 264mm - 278mm       | 252mm - 266mm        | 43mm           | 59mm       |
| 5    | 101mm       | 270mm - 284mm       | 258mm - 272mm        | 44mm           | 60mm       |
| 6    | 103mm       | 276mm - 290mm       | 264mm - 278mm        | 45mm           | 61mm       |
| 7    | 105mm       | 282mm - 296mm       | 270mm - 284mm        | 46mm           | 62mm       |
| 8    | 107mm       | 288mm - 302mm       | 276mm - 290mm        | 47mm           | 63mm       |
| 9    | 109mm       | 294mm - 308mm       | 282mm - 296mm        | 48mm           | 64mm       |

# MEN'S ACCOMODATION CHART

## MEN'S NARROW LOWLINE

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 5    | 85mm        | 224mm - 234mm       | 222mm - 232mm        | 32mm           | 58mm       |
| 6    | 88mm        | 230mm - 240mm       | 228mm - 238mm        | 33mm           | 59mm       |
| 7    | 89mm        | 236mm - 246mm       | 234mm - 244mm        | 34mm           | 60mm       |
| 8    | 91mm        | 242mm - 252mm       | 240mm - 250mm        | 35mm           | 61mm       |
| 9    | 93mm        | 249mm - 259mm       | 247mm - 257mm        | 36mm           | 62mm       |
| 10   | 95mm        | 255mm - 265mm       | 253mm - 263mm        | 37mm           | 63mm       |
| 11   | 97mm        | 263mm - 273mm       | 261mm - 271mm        | 38mm           | 64mm       |
| 12   | 99mm        | 269mm - 279mm       | 267mm - 277mm        | 39mm           | 65mm       |

## MEN'S NARROW

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 5    | 85mm        | 235mm - 248mm       | 228mm - 241mm        | 35mm           | 58mm       |
| 6    | 88mm        | 241mm - 254mm       | 235mm - 248mm        | 36mm           | 59mm       |
| 7    | 89mm        | 248mm - 260mm       | 241mm - 254mm        | 37mm           | 60mm       |
| 8    | 91mm        | 254mm - 266mm       | 248mm - 260mm        | 38mm           | 61mm       |
| 9    | 93mm        | 260mm - 273mm       | 254mm - 266mm        | 39mm           | 62mm       |
| 10   | 95mm        | 266mm - 279mm       | 260mm - 273mm        | 40mm           | 63mm       |
| 11   | 97mm        | 273mm - 287mm       | 266mm - 279mm        | 41mm           | 64mm       |
| 12   | 99mm        | 279mm - 293mm       | 273mm - 287mm        | 42mm           | 65mm       |

## MEN'S STANDARD LOWLINE

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 5    | 94mm        | 236mm - 246mm       | 234mm - 244mm        | 32mm           | 61mm       |
| 6    | 95mm        | 242mm - 252mm       | 240mm - 250mm        | 33mm           | 62mm       |
| 7    | 97mm        | 249mm - 259mm       | 246mm - 256mm        | 34mm           | 63mm       |
| 8    | 99mm        | 255mm - 265mm       | 253mm - 263mm        | 35mm           | 64mm       |
| 9    | 102mm       | 263mm - 273mm       | 261mm - 271mm        | 36mm           | 65mm       |
| 10   | 105mm       | 269mm - 279mm       | 267mm - 277mm        | 37mm           | 66mm       |
| 11   | 108mm       | 275mm - 285mm       | 273mm - 283mm        | 38mm           | 67mm       |
| 12   | 112mm       | 282mm - 292mm       | 280mm - 290mm        | 39mm           | 68mm       |

## MEN'S STANDARD

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 5    | 94mm        | 248mm - 260mm       | 241mm - 254mm        | 35mm           | 61mm       |
| 6    | 95mm        | 254mm - 266mm       | 248mm - 260mm        | 36mm           | 62mm       |
| 7    | 97mm        | 260mm - 273mm       | 254mm - 266mm        | 37mm           | 63mm       |
| 8    | 99mm        | 266mm - 279mm       | 260mm - 273mm        | 38mm           | 64mm       |
| 9    | 102mm       | 273mm - 287mm       | 266mm - 279mm        | 39mm           | 65mm       |
| 10   | 105mm       | 279mm - 293mm       | 273mm - 287mm        | 40mm           | 66mm       |
| 11   | 108mm       | 287mm - 299mm       | 279mm - 293mm        | 41mm           | 67mm       |
| 12   | 112mm       | 293mm - 306mm       | 287mm - 299mm        | 42mm           | 68mm       |

## MEN'S WIDE LOWLINE

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 5    | 99mm        | 246mm - 256mm       | 238mm - 248mm        | 32mm           | 61mm       |
| 6    | 102mm       | 253mm - 263mm       | 245mm - 255mm        | 33mm           | 62mm       |
| 7    | 105mm       | 259mm - 269mm       | 252mm - 262mm        | 34mm           | 63mm       |
| 8    | 108mm       | 265mm - 275mm       | 259mm - 269mm        | 35mm           | 64mm       |
| 9    | 112mm       | 272mm - 282mm       | 266mm - 276mm        | 36mm           | 65mm       |
| 10   | 115mm       | 278mm - 288mm       | 272mm - 282mm        | 37mm           | 66mm       |
| 11   | 118mm       | 286mm - 296mm       | 278mm - 288mm        | 38mm           | 67mm       |
| 12   | 120mm       | 294mm - 302mm       | 284mm - 294mm        | 39mm           | 68mm       |

## MEN'S WIDE

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 5    | 99mm        | 258mm - 270mm       | 248mm - 254mm        | 35mm           | 61mm       |
| 6    | 102mm       | 264mm - 277mm       | 252mm - 260mm        | 36mm           | 62mm       |
| 7    | 105mm       | 270mm - 283mm       | 254mm - 264mm        | 37mm           | 63mm       |
| 8    | 108mm       | 277mm - 290mm       | 264mm - 270mm        | 38mm           | 64mm       |
| 9    | 112mm       | 283mm - 296mm       | 270mm - 277mm        | 39mm           | 65mm       |
| 10   | 115mm       | 290mm - 302mm       | 277mm - 283mm        | 40mm           | 66mm       |
| 11   | 118mm       | 299mm - 311mm       | 285mm - 295mm        | 41mm           | 67mm       |
| 12   | 120mm       | 308mm - 320mm       | 295mm - 301mm        | 42mm           | 68mm       |

## MEN'S DEEP


| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 5    | 94mm        | 266mm - 279mm       | 260mm - 273mm        | 44mm           | 61mm       |
| 6    | 95mm        | 273mm - 287mm       | 266mm - 279mm        | 45mm           | 62mm       |
| 7    | 97mm        | 279mm - 293mm       | 273mm - 287mm        | 46mm           | 63mm       |
| 8    | 99mm        | 287mm - 299mm       | 279mm - 293mm        | 47mm           | 64mm       |
| 9    | 102mm       | 293mm - 306mm       | 287mm - 299mm        | 48mm           | 65mm       |
| 10   | 105mm       | 299mm - 312mm       | 293mm - 306mm        | 49mm           | 66mm       |
| 11   | 108mm       | 306mm - 318mm       | 299mm - 312mm        | 50mm           | 67mm       |
| 12   | 112mm       | 308mm - 320mm       | 306mm - 314mm        | 51mm           | 68mm       |

## MEN'S WIDE & DEEP

| SIZE | JOINT WIDTH | JOINT CIRCUMFERENCE | INSTEP CIRCUMFERENCE | TOE BOX DEPTHS | HEEL WIDTH |
|------|-------------|---------------------|----------------------|----------------|------------|
| 5    | 99mm        | 266mm - 280mm       | 259mm - 273mm        | 44mm           | 61mm       |
| 6    | 102mm       | 274mm - 288mm       | 265mm - 279mm        | 45mm           | 62mm       |
| 7    | 105mm       | 282mm - 296mm       | 272mm - 286mm        | 46mm           | 63mm       |
| 8    | 108mm       | 290mm - 304mm       | 279mm - 293mm        | 47mm           | 64mm       |
| 9    | 112mm       | 298mm - 312mm       | 285mm - 299mm        | 48mm           | 65mm       |
| 10   | 115mm       | 306mm - 320mm       | 292mm - 306mm        | 49mm           | 66mm       |
| 11   | 118mm       | 314mm - 328mm       | 298mm - 312mm        | 50mm           | 67mm       |
| 12   | 120mm       | 322mm - 336mm       | 304mm - 318mm        | 51mm           | 68mm       |

## LEATHER/SOFTEE LEATHER

Leather is durable with strength and flexibility. Softee leather is smooth and supple.



## NUBUCK

Buffed smooth nap full-grain leather.



## PATENT

Firm with high-gloss film finish.




## MILLED

Soft and supple with print finish.




## SUEDE

Go Rain water resistant, medium softness.




## BRUISER

Heavyweight supple with soft handle.




## NEOPRENE

Breathable stretch material.




## AIRNET

3d spacer fabric, excellent breathability and air circulation.



## HOUSE SHOE FELT

Comfortable, naturally insulating and breathable.



## SOLING OPTIONS



Pyramid



Nova



Oporto



Quinas



Ribbed PU



Commando



Ripple



TCR Oil Resistant

## SOLE PATTERN OPTIONS



ULLI EVA



WELLAN



RIBBED RUBBER



TYRE TRACK RUBBER

## LACING/UPPER ALTERNATIVES



Elastic Toggle



Twin Grip Velcro



Soft Top Collar

# TERMS AND CONDITIONS

## FOOTWEAR

All of our footwear can be made with a diabetic specification of soft leather, rim toe puff and Superplush antibacterial lining.

—

All of our styles can be made on any of our 3 lasts. Your order will be made on the last pictured in the catalogue unless alternative is specified.

—

EVA is used for Standard sole and heels and through wedge soles and soles requiring raises, sockets, etc.

—

Hidden depth sole unit is used on Nova last styles unless alternative requested with order or width additions exceed 6mm or floats, sockets are required.

—

Sole used will be as catalogue image unless an alternative is specified.

—

First trial fitting is **free of charge**, subsequent fittings may be charged for. Accurate measurements are essential when placing your order as excessive adaption following trial fit which require new uppers, etc. will be charged for.

—

Documentation relating to the confirmation of orders should be clearly marked '**confirmation**' to avoid duplication. Duplicated orders through originator's error will be chargeable.

—

The company reserves the right to alter specification or withdraw styles without prior notice. Specification changes are only undertaken when product improvements are identified or a component material is no longer available.

—

Colour reproduction - We make every effort to ensure the colour you see in our catalogue is accurate, however, colours can be different in photographs and shades can vary over the life of our catalogue.

## READY-MADE/STOCK ITEMS

All of our ready-made footwear - Podartis, Podowell, and Lymph are returnable for a full refund; they must be unworn and in their original undamaged packaging.

## DESPATCH

UK despatch is **next day service**. Next day timed delivery is available but will incur extra cost.

—

An invoice will be raised for all orders not returned for completion within 90 days.

—

All goods remain the property of Reed Medical Ltd until paid for in full.

## RETURNS

All our Reeco footwear is made to order therefore non-returnable with the exception only of a manufacturing fault of our making.

## CATALOGUES/WEBSITE

All of our catalogues and brochures are available to view or download from our website: [www.reedmedical.co.uk](http://www.reedmedical.co.uk).

Our website is a live document and we constantly feed it with information regarding all of our new styles, current news items, etc. Please visit it regularly or follow us on LinkedIn or Facebook where you will receive updates/information regularly.

—

Should you require further information or wish to speak with our Customer Care and Technical team please contact us on **01254 503333**.

—

**N.B. Please note that we do not employ an Orthotist or offer medical advice.**



See our full range of Kids' footwear by  
downloading the catalogue from our website  
[www.reedmedical.co.uk](http://www.reedmedical.co.uk)

Scan me!



# reed

bespoke footwear

[www.reedmedical.co.uk](http://www.reedmedical.co.uk)

Tel. 01254 503333

[admin@reedmedical.co.uk](mailto:admin@reedmedical.co.uk)

Scan me!



FOLLOW US ON

