

PODARTIS FOOTWEAR SIZING GUIDE

CHART A												
SIZES	36	37	38	39	40	41	42	43	44	45	46	47
JOINT WIDTH	90	92	94	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	243	248	253	258	263	268	273	278	283	288	293	298
TOE BOX HEIGHT	32	33	34	35	36	37	38	39	40	41	42	43

PRODUCT RANGE
DEAMBULO L
DEAMBULO X

CHART B							
SIZES	35	36	37	38	39	40	41
JOINT WIDTH	87	89	90	92	93	95	96
JOINT CIRCUMFERENCE	238	241	246	253	258	261	268
TOE BOX HEIGHT	31	32	33	34	35	36	37

PRODUCT RANGE
FLY
DEAMBULO SPORT

CHART C							
SIZES	35	36	37	38	39	40	41
JOINT WIDTH	87	89	90	92	93	95	96
JOINT CIRCUMFERENCE	241	246	251	256	261	266	271
TOE BOX HEIGHT	36	37	38	39	40	41	42

PRODUCT RANGE
FANCY
ITALIA
EVO MIDNIGHT

CHART D							
SIZES	35	36	37	38	39	40	41
JOINT WIDTH	85	87	88	90	91	93	94
JOINT CIRCUMFERENCE	241	246	251	256	261	266	271
TOE BOX HEIGHT	33	34	35	36	37	38	39

PRODUCT RANGE
ADELLE VELCRO

PODARTIS FOOTWEAR SIZING GUIDE

CHART E									
SIZES	39	40	41	42	43	44	45	46	47
JOINT WIDTH	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	261	266	271	276	281	286	291	296	301
TOE BOX HEIGHT	34	34	34	35	36	37	38	39	40

PRODUCT RANGE
HERO EVO
ONE EVO MAN

CHART F													
SIZES	35	36	37	38	39	40	41	42	43	44	45	46	47
JOINT WIDTH	87	89	91	93	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	241	246	251	256	261	266	271	276	281	286	291	296	301
TOE BOX HEIGHT	34	34	34	34	34	34	34	35	36	37	38	39	40

PRODUCT RANGE
FREEDOM EVO

CHART G									
SIZES	39	40	41	42	43	44	45	46	47
JOINT WIDTH	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	261	266	271	276	281	286	291	295	301
TOE BOX HEIGHT	36	37	38	39	40	41	42	43	44

PRODUCT RANGE
DCS SPORT
X-DIAB 14

CHART H									
SIZES	39	40	41	42	43	44	45	46	47
JOINT WIDTH	96	97	99	100	102	103	105	106	108
JOINT CIRCUMFERENCE	256	261	266	271	276	281	286	291	296
TOE BOX HEIGHT	34	35	36	37	38	39	40	41	42

PRODUCT RANGE
GIOTTO
BOTERO 11 VELCRO

PODARTIS FOOTWEAR SIZING GUIDE

CHART I								
SIZES	35	36	37	38	39	40	41	42
JOINT WIDTH	87	89	90	92	93	95	96	97
JOINT CIRCUMFERENCE	251	256	261	266	271	276	281	286
TOE BOX HEIGHT	37	38	39	40	41	42	43	44

PRODUCT RANGE
LADY ALFA

CHART J								
SIZES	39	40	41	42	43	44	45	46
JOINT WIDTH	96	97	99	100	102	103	105	106
JOINT CIRCUMFERENCE	271	276	281	286	291	296	301	306
TOE BOX HEIGHT	39	40	41	42	43	44	45	46

PRODUCT RANGE
GIOTTO VELCRO GV

CHART K											
SIZES	36	37	38	39	40	41	42	43	44	45	46
JOINT WIDTH	93	95	96	98	99	101	102	104	105	107	108
JOINT CIRCUMFERENCE	262	267	272	277	282	287	292	297	302	307	312
TOE BOX HEIGHT	39	40	41	42	43	44	45	46	47	48	49

PRODUCT RANGE
X-DIAB 16