NuStride

Kids Collection

2026



Page 4 Dok



Page 11 Jamie



Page 13 Taylor

MADE IN BRITAIN

reed by eqwal.



Pictured: Rose

WELCOME

to the 9th edition of the Reed Children's Footwear Catalogue

Based in the North West of England, our experienced technicians have been manufacturing specialist children's footwear for over 60 years. By integrating traditional craftsmanship with modern manufacturing processes and technology, we produce footwear designed to support clinical requirements while maintaining a high standard of durability and appearance.

CONTENTS

Kidztuff Collection 21-40	4-5
School Collection 21-40	6-7
Leisure Collection 21-40	8-12
Insoles	12
Tiny Toes 15-20	13
Anti-Toe Walker Collection	14
Embroidery Options	15
Colour / Material Chart	16
Options	17
Accommodation Chart	18
Order Form	19

reed by eqwal.

MADE IN BRITAIN

nqa ISO 13485 MEDICAL DEVICES UKAS MANAGEMENT SYSTEMS 0015

214424

nqa ISO 9001 QUALITY MANAGEMENT UKAS MANAGEMENT SYSTEMS 0015

214423

KIDZTUFF® 21-40

All the footwear on these two pages are made to our Kidztuff® specification. Shown on our natural last.



SATURN
Shown in: Black bruiser, Nubuck collar, commando sole.



MERCURY
Shown in: Black bruiser, commando sole.



PLUTO
Shown in: Black bruiser anti-scuff mudguard & counter, commando sole.



DOK
Shown in: Black bruiser, commando/yellow stitch.

Available with a medial zip option. Note: Medially extended stiffeners are not available with zip option.

- Strong supple upper leather
- Full leather linings
- Double strength medially extended heel stiffeners
- Foam padded quarters
- Commando or EVA/rubber outsole



ROBIN
Shown in: Black bruiser, through wedge sole.



ROSE
Shown in: Black bruiser, black embroidery, commando sole/yellow stitch.

Available with a medial zip option.
Note: Medially extended stiffeners are not available with zip option.



STAR
Shown in: Black bruiser, black daisy, commando sole.



IRIS
Shown in: Black bruiser, patent heart, through wedge sole.

SCHOOL COLLECTION 21-40

All the footwear on these two pages are made to our standard specification. Shown on our classic last.



RAVEN
Shown in: Black calf, Nubuck collar, sole & heel.



DODGER
Shown in: Black calf, Nubuck collar, sole & heel.



GRACE
Shown in: Black calf, black embroidery, commando sole, yellow stitch

Available with a medial zip option.
Note: Medially extended stiffeners are not available with zip option.



DARCEY
Shown in: Black calf, red patent heart, through wedge sole.

- Calf leathers
- Choice of leather or superplush lining
- Single strength standard stiffeners
- Commando or EVA/rubber outsole
- EVA sole/heel unit



JAY
Shown in: Black calf with medial zip, commando sole/yellow stitch.

Note: Medially extended stiffeners are not available with zip option.



BLOSSOM
Shown in: Black calf, daisy quarters, commando sole.



SAPPHIRE
Shown in: Black calf, sole & heel.



BELL
Shown in: Black calf, black pebble, through wedge sole.

LEISURE 21-40

All the footwear on these pages are made to our standard specification.
Shown on our classic last.



FREYA

Shown in: Black calf, red/green embroidery, commando sole.

Available with a medial zip option.
Note: Medially extended stiffeners are not available with zip option.



COUGAR

Shown in: Honey Nubuck, brown leather collar, commando sole.



BLOSSOM

Shown in: Purple calf with purple animal print suede quarters, commando sole.



JAY

Shown in: Brown softee with medial zip, commando sole.

Note: Medially extended stiffeners are not available with zip option.

- Choice of material and colour
- Leather or absorber lining
- Single strength stiffeners
- Commando or EVA/rubber outsole



TAYLOR
Shown in: Lilac animal print suede, through wedge sole.



TIGER
Shown in: White softee, through wedge sole.



SPYDER
Shown in: White milled, pink animal print flash, through wedge sole.



PACER
Shown in: Grey milled, white leather flash, through wedge sole.

LEISURE 21-40

All the footwear on these pages are made to our standard specification. Shown on our classic last.



VYPER
Shown in: Navy milled, tan Nubuck flash, through wedge sole.



REESE
Shown in: Red suede, red-black airnet, through wedge sole.



TREK
Shown in: Red softee, white contrast stitch, through wedge sole.



DUNE
Shown in: Tan softee (Velcro only), EVA sole/heel.

- Choice of material and colour
- Leather or absorber lining
- Single strength stiffeners
- Commando or EVA/rubber outsole



JAMIE
Shown in: Dark green suede, white leather flash, through wedge sole.



JORDAN
Shown in: Black & white softee, through wedge sole.



HARLEY
Shown in: Black softee, black suede, white leather flash, through wedge sole.



PARKER
Shown in: Royal blue airnet, blue suede, reflective flash, through wedge sole.

LEISURE 21-40

All the footwear on these pages are made to our standard specification. Shown on our classic last.



OLLIE
Shown in: White softee, navy suede, red flash, through wedge sole.



BELL
Shown in: White softee, mint milled, through wedge sole.

INSOLES

Choose one of our insoles to be added to your footwear order (optional, added cost).

ANTI-PRONATION INSOLE

The anti-pronation foot orthosis is designed to resist excessive pronation.

It incorporates a medium-density base, a microfibre top cover, and a 3° intrinsic medial heel post.



NEURO INSOLE

The neurological foot orthosis incorporates a deep heel cup, an enhanced transverse arch, and a second-to-fifth toe crest.



TINY TOES COLLECTION 15-20

- Extended opening
- Single strength stiffener
- Choice of leather
- Through wedge sole



RAVEN
Shown in: Red calf, Nubuck collar, through wedge sole.



SQUEAK
Shown in: Natural milled, cream suede collar, through wedge sole. Teddy embroidered tongue.



TAYLOR
Shown in: Soft pink saffiano, through wedge sole.



VYPER
Shown in: Navy milled, white milled flash through wedge sole.

Embroidery not available with this style.

ANTI-TOE WALKER COLLECTION 21-40

All the boots on this page are made to our Kidztuff specification with the addition of a padded tongue, a 30mm increase in back height and carbon fibre insert, made on our Natural Toe last.



REGGAE
Shown in: Black bruiser, Nubuck collar, through wedge sole.



FOX
Shown in: Black bruiser, commando sole.

Available with a medial zip option.
Note: Medially extended stiffeners are not available with zip option.



KEPLER
Shown in: Black bruiser, anti-scuff mudguard & counter, commando sole.



RORY
Shown in: Black bruiser, through wedge sole.

EMBROIDERY OPTIONS

Choose an embroidery design to be placed on the tongue or outside quarter of your boot (optional extra).



TEDDY
Available 15 - 40



UNICORN
Available 15 - 40



HEART
Available 15 - 40



FOOTBALL
Available 15 - 40



SNAKE
Available 15 - 40



BUTTERFLY
Available 18 - 40



DINOSAUR
Available 18 - 40



CONTROLLER
Available 18 - 40



BALLERINA
Available 21 - 40



SINGLE ROSE
Available 15 - 40



DOUBLE ROSE
Available 21 - 40



ROSES
Available 31 - 40

COLOUR CHART

We offer a comprehensive selection of upper materials, including Bruiser leather, Calf/Softee leather, Nubuck, Suede, Airnet, Anti-Scuff and Decorative finishes. These materials are selected for their performance characteristics, including durability, breathability, and adaptability.

LEATHER / SOFTEE LEATHER

Leather is durable with strength & flexibility



BRUISER LEATHER

Heavyweight leather supple with soft handle



SUEDE

Go rain water resistant, medium softness



NUBUCK

Buffed smooth nap full grain leather



PATENT

Firm with high gloss film finish



ANIMAL PRINT

Medium softness suede, gloss film print



DAISY



AIRNET

3D spacer fabric, excellent breathability & air circulation




MILLED

Soft and supple with print finish



SAFFIANO

A full grain, with a glossy finish and cross-hatch pattern



SPARKLE

Soft suede with metallic glitter



OPTIONS

Our range is built around four anatomically considered last shapes — Classic, Natural/Straight, Reverse, and AFO — developed to accommodate a variety of foot presentations and orthotic needs. Styles are available in sizes 15–40 and across five width fittings to support an accurate and consistent fit.

LASTS:

Children's lasts include a 3mm insole allowance



Classic/Round



AFO



Natural/Straight



Reverse

STIFFENERS:

Available in Single or Double Strength



Standard



Medially Extended




High & Extended



High


SOLES:



Commando

(Hardwearing Rubber)


Black or Caramel



Sole & Heel

(Hardwearing EVA)

Black only



Through Wedge

(EVA Wedge / Rubber Outsole)

Black or White

Kids Last Accommodation Chart

Kids Extra Extra Narrow

Size	Joint Width	Joint Circ.	Instep Circ.
15	43mm	115mm	125mm
16	44mm	120mm	130mm
17	46mm	125mm	135mm
18	47mm	130mm	140mm
19	49mm	135mm	145mm
20	52mm	140mm	150mm
21	53mm	145mm	155mm
22	55mm	150mm	160mm
23	56mm	155mm	165mm
24	58mm	160mm	170mm
25	59mm	165mm	175mm
26	60mm	170mm	180mm
27	62mm	175mm	185mm
28	63mm	180mm	190mm
29	65mm	185mm	195mm
30	67mm	190mm	200mm
31	69mm	195mm	205mm
32	70mm	200mm	210mm
33	71mm	205mm	215mm
34	73mm	210mm	220mm
35	75mm	215mm	225mm
36	76mm	220mm	230mm
37	78mm	225mm	235mm
38	80mm	230mm	240mm
39	82mm	235mm	245mm
40	84mm	240mm	250mm

Kids Extra Narrow

Size	Joint Width	Joint Circ.	Instep Circ.
15	47mm	120mm	130mm
16	48mm	125mm	135mm
17	50mm	130mm	140mm
18	51mm	135mm	145mm
19	53mm	140mm	150mm
20	56mm	145mm	155mm
21	57mm	150mm	160mm
22	59mm	155mm	165mm
23	60mm	160mm	170mm
24	62mm	165mm	175mm
25	63mm	170mm	180mm
26	64mm	175mm	185mm
27	66mm	180mm	190mm
28	67mm	185mm	195mm
29	69mm	190mm	200mm
30	71mm	195mm	205mm
31	73mm	200mm	210mm
32	74mm	205mm	215mm
33	75mm	210mm	220mm
34	77mm	215mm	225mm
35	79mm	220mm	230mm
36	80mm	225mm	235mm
37	82mm	230mm	240mm
38	84mm	235mm	245mm
39	86mm	240mm	250mm
40	88mm	245mm	255mm

Kids Narrow

Size	Joint Width	Joint Circ.	Instep Circ.
15	51mm	130mm	140mm
16	52mm	135mm	145mm
17	54mm	140mm	150mm
18	55mm	145mm	155mm
19	58mm	150mm	160mm
20	60mm	155mm	165mm
21	61mm	160mm	170mm
22	63mm	165mm	175mm
23	64mm	170mm	180mm
24	66mm	175mm	185mm
25	67mm	180mm	190mm
26	68mm	185mm	195mm
27	70mm	190mm	200mm
28	71mm	195mm	205mm
29	73mm	200mm	210mm
30	75mm	205mm	215mm
31	77mm	210mm	220mm
32	78mm	215mm	225mm
33	79mm	220mm	230mm
34	81mm	225mm	235mm
35	83mm	230mm	240mm
36	84mm	235mm	245mm
37	86mm	240mm	250mm
38	88mm	245mm	255mm
39	90mm	250mm	260mm
40	92mm	255mm	265mm

Kids Standard


Size	Joint Width	Joint Circ.	Instep Circ.
15	54mm	135mm	145mm
16	56mm	140mm	150mm
17	57mm	145mm	155mm
18	59mm	150mm	160mm
19	60mm	155mm	165mm
20	62mm	160mm	170mm
21	63mm	165mm	175mm
22	65mm	170mm	180mm
23	66mm	175mm	185mm
24	68mm	180mm	190mm
25	69mm	185mm	195mm
26	71mm	190mm	200mm
27	72mm	195mm	205mm
28	74mm	200mm	210mm
29	75mm	205mm	215mm
30	77mm	210mm	220mm
31	78mm	215mm	225mm
32	80mm	220mm	230mm
33	81mm	225mm	235mm
34	83mm	230mm	240mm
35	84mm	235mm	245mm
36	87mm	240mm	250mm
37	89mm	245mm	255mm
38	91mm	250mm	260mm
39	92mm	255mm	265mm
40	94mm	260mm	270mm

Kids Wide

Size	Joint Width	Joint Circ.	Instep Circ.
15	56mm	140mm	150mm
16	58mm	145mm	155mm
17	59mm	150mm	160mm
18	61mm	155mm	165mm
19	63mm	160mm	170mm
20	65mm	165mm	175mm
21	67mm	170mm	180mm
22	69mm	175mm	185mm
23	71mm	180mm	190mm
24	73mm	185mm	195mm
25	74mm	190mm	200mm
26	75mm	195mm	205mm
27	76mm	200mm	210mm
28	78mm	205mm	215mm
29	79mm	210mm	220mm
30	81mm	215mm	225mm
31	83mm	220mm	230mm
32	85mm	225mm	235mm
33	87mm	230mm	240mm
34	89mm	235mm	245mm
35	90mm	240mm	250mm
36	92mm	245mm	255mm
37	93mm	250mm	260mm
38	95mm	255mm	265mm
39	96mm	260mm	270mm
40	98mm	265mm	275mm

Kids AFO

Size	Joint Width	Joint Circ.	Instep Circ.
15	56mm	145mm	160mm
16	58mm	150mm	165mm
17	59mm	155mm	170mm
18	61mm	160mm	175mm
19	63mm	165mm	180mm
20	65mm	170mm	185mm
21	67mm	175mm	190mm
22	69mm	180mm	195mm
23	71mm	185mm	200mm
24	73mm	190mm	205mm
25	74mm	195mm	210mm
26	75mm	200mm	215mm
27	76mm	205mm	220mm
28	78mm	210mm	225mm
29	79mm	215mm	230mm
30	81mm	220mm	235mm
31	83mm	225mm	240mm
32	85mm	230mm	245mm
33	87mm	235mm	250mm
34	89mm	240mm	255mm
35	90mm	245mm	260mm
36	92mm	250mm	265mm
37	93mm	255mm	270mm
38	95mm	260mm	275mm
39	96mm	265mm	280mm
40	98mm	270mm	285mm



An AFO last is a standard last modified to accommodate a patient wearing an AFO. Adjustments include an additional 10 mm in length, 15 mm in instep circumference, 8 mm in seat depth and 10 mm in seat width.

Paediatric Footwear Order Form

reed by eqwal.

Order Number:_____ Date:_____

Patient Name:_____

Clinician Name:_____

Footwear Options

Style Name:_____ Trim Colour:_____

Colour/Material:_____ Lining:_____

☐ Lace ☐ Velcro

Boot Height: ☐ Standard ☐ Bespoke _____mm Top Circ (A measure) L_____mm R_____mm

Last Options

Size: Left_____ Right_____

☐ Round ☐ Natural ☐ Afo ☐ Reverse

LEFT	RIGHT
☐ X X Narrow	☐ X X Narrow
☐ X Narrow	☐ X Narrow
☐ Narrow	☐ Narrow
☐ Standard	☐ Standard
☐ Wide	☐ Wide

Sole

☐ Commando ☐ Through Wedge

☐ Sole & Heel

Additional Sole Instructions:

_____

_____

_____

_____

_____

Stiffeners

☐ Standard Strength ☐ Double Strength

☐ Padded

Insoles

☐ PL930 Disport ☐ Anti-Pronation

☐ PL910 Star ☐ Neurological

Trial Fit? ☐

Last Adaptations

LEFT:

_____

_____

_____

RIGHT:

_____

_____

_____

Additional information:

+44 01254 503333 www.reedmedical.co.uk orthotics@reedmedical.co.uk

reed
by eqwal.

NuStride
Kids Collection

www.reedmedical.co.uk

Email: orthotics@reedmedical.co.uk

Tel: 01254 503333

Greenbank Mill, Gladstone St,
Blackburn, Lancashire, BB1 3ES

FOLLOW US ON

Scan to view our brochures online!

Kids' Catalogue

Adults' Catalogue

Off The Shelf Catalogue